

BUILT IN SPORT. MADE FOR LIFE.™

PRODUCT SHOWCASE



MUELLER



BUILT IN SPORT. MADE FOR LIFE.®



For over 60 years, Mueller Sports Medicine has stayed close to the changing needs of athletes, be they weekend warriors or international superstars.

The mission never changes:

Help athletes stay in the game and help them perform at their peak.

We Are Sports Medicine®



WE ARE SPORTS MEDICINE

TABLE OF CONTENTS

5-11 TEAM TAPES – *BULK TAPES*

Rigid Adhesive Tapes
Elastic Adhesive Tapes
Cohesive Tapes
Specialty Tapes
Underwrap & Wraps

12 TAPING ACCESSORIES

13-17 ATHLETIC TRAINER KITS

18-21 RECOVERYCARE

REVIVE
Recovery Tubs
Therm X

22-24 COLD/HOT THERAPY

25-28 SPORT ACCESSORIES

29-33 ATHLETE PROTECTION

34-43 ANKLE BRACE & SUPPORTS

Advanced Ankle Supports
Wraparound Ankle Supports
Ankle Supports
Rigid Ankle Braces
Knee Straps

44-57 KNEE BRACE & SUPPORTS

Hinged Knee Braces
Knee Stabilizers & Braces
Knee Supports
Knee Sleeves

58-65 ELBOW BRACE & SUPPORTS

Adjustable Elbow Supports
Elbow Supports
Tennis Elbow Supports

66-70 WRIST BRACE & SUPPORTS

Wrist Braces & Stabilizers
Wrist Supports

71-73 BACK BRACE & SUPPORTS

74-77 GRADUATED COMPRESSION

78-80 SPECIALTY BRACES & SUPPORTS

81-82 QUENCH® GUM

RIGID ADHESIVE TAPE - *BULK TAPE* SOLD IN CASES

Traditional approach to providing support to joints and body parts, be it for injury prevention or injury management. Used when there is a need to restrict the range of motion of a given joint. Proper taping techniques can aid in the treatment of a wide range of painful musculoskeletal conditions such as lateral epicondylitis (tennis elbow), plantar fasciitis, osteoarthritis, anterior knee pain, patellofemoral pain syndrome, wry neck, shoulder subluxation and shoulder impingement syndrome. Most common is a rigid ankle taping technique that helps improve joint position sense.


Adhesive tapes can be manufactured in two main way, porous and non-porous. Porous tape is manufactured in a manner which facilitates breathability without compromising strength. Sweat vapor moves through the tape which provides greater comfort, especially to those with sensitive skin. Non-porous tape is the traditional approach to tape, where the adhesive coats the entirety of the tape surface. Both processes are beneficial based on your preference.

Porous

	Perform™ Plus Tape	Perform™ Tape	Athletic Care® Porous Tape
Size/Case	1.5" x 15 yds / 32 rolls 2" x 15 yds / 24 rolls	1.5" x 15 yds / 32 rolls 2" x 15 yds / 24 rolls	1.5" x 15 yds / 32 rolls
Fabric	100% Cotton, Woven	100% Cotton, Woven	100% Cotton, Zinc oxide
Water Resistant	✘	✘	✘
Not Made With Natural Rubber Latex			✘
Edge Style	Serrated	Straight	Straight

Non-Porous

	MTape®	MTape® Team Colors	Athletic Care® Athletic Trainers' Tape
Size/Case	1" x 10 yds / 48 rolls 1.5" x 15 yds / 32 rolls 2" x 15 yds / 24 rolls	1.5" x 10 yds / 32 rolls	1.5" x 15 yds / 32 rolls
Fabric	100% Cotton, Zinc oxide	100% Cotton, Woven, Zinc Oxide	100% Cotton, Zinc oxide
Water Resistant	✘	✘	✘
Not Made With Natural Rubber Latex			✘
Edge Style	Straight	Straight	Straight



Colors Available



MTape®

Among the most popular tapes in the world, MTAPE® is Mueller's flagship tape.

- | | | |
|--|--|---|
| 130104 White 1" x 10 yds (48/cs)
MSRP: \$77.99 USD | 130105 White 1.5" x 15 yds (32/cs) | 130106 White 2" x 15 yds (24/cs) |
| 130824 Black 1.5" x 10 yds (32/cs)
MSRP: \$85.99 USD | 130821 Forest Green 1.5" x 10 yds (32/cs) | 130822 Scarlet 1.5" x 10 yds (32/cs) |
| 130820 Royal Blue 1.5" x 10 yds (32/cs) | 130826 Purple 1.5" x 10 yds (32/cs) | 130823 Gold 1.5" x 10 yds (32/cs) |
| 130825 Orange 1.5" x 10 yds (32/cs) | 130830 Pink 1.5" x 10 yds (32/cs) | 130828 Maroon 1.5" x 10 yds (32/cs) |
| 130829 Gray 1.5" x 10 yds (32/cs)
MSRP: \$85.99 USD | | |



Perform Plus™ Tape


Premium porous quality zinc oxide based tape for less skin irritation and increased breathability. Includes a higher thread count for increased strength and serrated edges for easy to tear.

- | | |
|--|---|
| 130182 White 1.5" x 15 yds (32/cs)
MSRP: \$85.99 USD | 130183 White 2" x 15 yds (24/cs) |
|--|---|

Perform™ Tape

A porous zinc oxide based tape for less skin irritation and increased breathability.

- | | |
|---|--|
| 23162 White 1.5" x 15 yds (32/cs)
MSRP: \$76.99 USD | 23163 White 2" x 15 yds (24/cs) |
|---|--|



Athletic Care® Athletic Trainers' Tape

Lightweight non-porous adhesive tape that is economically priced.

- | |
|--|
| 130888 White 1.5" x 15 yds (32/cs)
MSRP: \$56.99 USD |
|--|

Athletic Care® Porous Tape

Lightweight porous adhesive tape that is economically priced.

- | |
|--|
| 130399 White 1.5" x 15 yds 32/cs
MSRP: \$58.99 |
|--|

ELASTIC ADHESIVE TAPE - *BULK TAPE* SOLD IN CASES

Used when there is a need to help provide adjustable compression and support in the treatment of musculoskeletal injuries such as sprains, strains, dislocations, fractions and conditions such as arthritis. Ideal for holding items such as cold packs, padding or splints in place

	Stretch MTape® Premium <i>Heavyweight</i>	FlexLight™ Tape	Tear-Light™ Tape	Tear-Light™ Tape <i>Team Colors</i>
Size/Case	1" x 5 yds / 32 rolls 2" x 5 yds / 24 rolls 3" x 5 yds / 16 rolls	2" x 7.5 yds / 24 rolls 3" x 7.5 yds / 16 rolls	2" x 5 yds / 24 rolls 3" x 5 yds / 16 rolls 1.5" x 7.5 yds / 32 rolls 2" x 7.5 yds / 24 rolls 3" x 5 yds / 24 rolls	2" x 7.5 yds / 24 rolls
Fabric	100% Cotton	Cotton, Elastane	Cotton, Lycra Stretch, Woven	Cotton, Lycra Stretch, Woven
Water Resistant	x	x		
Not Made With Natural Rubber Latex		x	x	x
Edge Style	Straight	Straight	Straight	Straight



Stretch MTape® Premium


Maximum-strength heavyweight elastic adhesive tape.

21011 White 1" x 5 yds (32/cs)

21012 White 2" x 5 yds (24/cs)

21013 White 3" x 5 yds (16/cs)

MSRP: \$112.99 USD



FlexLight™ Tape

Lightweight elastic adhesive tape, thin and low profile for better conformity.

26438 White 2" x 7.5 yds (24/cs)

26448 White 3" x 7.5 yds (16/cs)

MSRP: \$74.99 USD

26447 Black 2" x 7.5 yds (24/cs)

26457 Black 3" x 7.5 yds (16/cs)

MSRP: \$79.99 USD



Tear-Light Tape®

Easy-to-tear elastic adhesive that provides superior support.

130622 White 2" x 5 yds (24/cs)

130623 White 3" x 5 yds (16/cs)

MSRP: \$60.99 USD

130626 White 1.5" x 7.5 yds (32/cs)

130632 White 2" x 7.5 yds (24/cs)

130633 White 3" x 7.5 yds (16/cs)

MSRP: \$71.99 USD

130642 Black 2" x 7.5 yds (24/cs)

130643 Black 3" x 7.5 yds (16/cs)

130681 Red 2" x 7.5 yds (24/cs)

130682 Pink 2" x 7.5 yds (24/cs)

130683 Green 2" x 7.5 yds (24/cs)

130684 Blue 2" x 7.5 yds (24/cs)

130685 Gold 2" x 7.5 yds (24/cs)

MSRP: \$76.99 USD



COHESIVE TAPE - *BULK TAPE* SOLD IN CASES

Stays in place while only sticking to itself, not to the skin while providing support and compression to body parts. Can be used to secure pads, bandages, cold packs and more. There are rigid cohesives that are used in similar applications as rigid adhesives, and elastic cohesive used in similar applications as adhesive elastics.

	Thor™ Tape	Recoil™ Tape	Tapewrap® Premium	Tapewrap® Premium Team Colors	MLastic®Tape
Size/Case	1.5" x 15 yds / 32 rolls	1" x 4 yds / 10 rolls 2" x 4 yds / 10 rolls 3.5" x 4 yds / 10 rolls	1" x 6 yds / 48 rolls 2" x 6 yds / 24 rolls 3" x 6 yds / 16 rolls	2" x 6 yds / 24 rolls 3" x 6 yds / 16 rolls	2" x 5 yds / 24 rolls 3" x 5 yds / 16 rolls
Fabric	Polyester, Cellulose	Thermoplastic Polyurethane	Viscose, Polyamide/Elastane	Cotton/Elastane	100% Cotton
Water Resistant	x	x			x
Not Made With Natural Rubber Latex					x
Edge Style	Straight	Straight	Straight	Straight	Straight



Thor™ Tape

Rigid cohesive tape that works with Mueller TapeWrap to provide maximum support.

26122 White 1.5" x 15 yds (32/cs)
MSRP: \$49.99 USD



Recoil™ Tape

Highly elastic, high-recovery, super-thin and water-resistant cohesive tape.

29010 White 1" x 4.4 yds (10/cs) 29310 Black 1" x 4.4 yds (10/cs)
MSRP: \$36.99 USD
29020 White 2" x 4.5 yds (10/cs) 29320 Black 2" x 4.4 yds (10/cs)
MSRP: \$61.99 USD
29035 White 3.5" x 4.4 yds (10/cs)
MSRP: \$99.99 USD



Tapewrap® Premium

Cohesive, lightweight and easy to tear tape. Also works as an underlayment for Thor™ Tape.

22057B White 1" x 6 yds 48/cs
MSRP: \$67.99 USD

24058B White 2" x 6 yds 24/cs 26059B White 3" x 6 yds 16/cs 24258B Black 2" x 6 yds 24/cs 26259B Black 3" x 6 yds 16/cs

MSRP: \$57.99 USD

24158 Grey 2" x 6 yds 24/cs 24160 Orange 2" x 6 yds 24/cs 24358 Beige 2" x 6 yds 24/cs 24458B Blue 2" x 6 yds 24/cs

24558B Red 2" x 6 yds 24/cs 24658B Pink 2" x 6 yds 24/cs 24758B Gold 2" x 6 yds 24/cs 24858B Green 2" x 6 yds 24/cs

24958B Purple 2" x 6 yds 24/cs
MSRP: \$58.99 USD



MLastic® Tape

Lightweight cohesive stretch tape.

130612 Beige 2" x 5 yds 24/cs
MSRP: \$80.99 USD
130613 Beige 3" x 5 yds 16/cs
MSRP: \$81.99 USD

KINESIOLOGY TAPE - BULK TAPE

	Kinesiology Tape Continuous Rolls	ProStrips®	PFTape®	PTape®
Size/Case	2" x 16.4' / 6 rolls 2" x 98.4' / 1 roll	2" x 5" / 24 strips (logo) 4" x 8" / 24 strips (logo) 2" x 10" / 1 roll (no logo) 4" x 10" / 1 roll (no logo) 6" x 10" / 1 roll (no logo)	REG / 56 apps	1.5" x 15 yds / 30 rolls
Fabric	Cotton, Polyurethane, Woven	Polyester Spunlace	Woven Rayon	Rayon
Water Resistant	✗		✗	✗
Not Made With Natural Rubber Latex	✗	✗	✗	✗
Edge Style	Straight	Straight	Straight	Serrated



Kinesiology Tape Continuous Rolls

Helps restore function to injured, sore or weak muscles and joints. Features a revolutionary wave pattern adhesive that moves with your skin and muscles as you heal.

27467 Beige 2" x 16.4' Rolls 6/tray

28147 Black 2" x 16.4' Rolls 6/tray

28277 Pink 2" x 16.4' Rolls 6/tray

27367 Blue 2" x 16.4' Rolls 6/tray

MSRP: \$65.99 US

27634 Beige 2" x 98.4'

27631 Black 2" x 98.4'

27633 Pink 2" x 98.4'

27632 Blue 2" x 98.4'

MSRP: \$39.99 USD

SPECIALTY TAPES - *BULK TAPE* SOLD IN CASES



ProStrips®

Easy to use pressure-sensitive adhesive mesh that helps protect against abrasions, turf burns, blisters and calluses. Great for racquet sports. Works well with PTape®

No logo

010972 2" x 10 yd Roll
MSRP: \$27.99 USD

010974 4" x 10 yd Roll
MSRP: \$38.99 USD

010976 6" x 10 yd Roll
MSRP: \$50.99 USD

Logo

010952 2" x 5" Pre-Cut
MSRP: \$23.99 USD

010954 4" x 8" Pre-Cut
MSRP: \$33.99 USD

PFTape®

PFTape® System for Plantar Fasciitis is a patented Mueller exclusive rigid tape used for treating heel and arch pain.

25337 Beige REG 56/box
MSRP: \$108.99 USD

PTape®

Rigid strapping tape that works great with ProStrips® for patellofemoral/McConnell taping methods.

21546 Beige 1.5" x 15 yds 30/cs (5 inner boxes of 6)
MSRP: \$306.99 USD

UNDERWRAP & WRAPS - *BULK* SOLD IN CASES



MWrap®

Pre-taping foam underwrap helps protect the skin from tape.

130701 Natural 2.75" x 30 yd Roll 12/cs
MSRP: \$44.99 USD

130704 Big Red 2.75" x 30 yd Roll 48/cs

130706 Big Gold 2.75" x 30 yd Roll 48/cs

130708 Big Blue 2.75" x 30 yd Roll 48/cs

130712 Big Purple 2.75" x 30 yd Roll 48/cs

130714 Big Pink 2.75" x 30 yd Roll 48/cs

130719 Big Aqua 2.75" x 30 yd Roll 48/cs

130722 Big Sunburst 2.75" x 30 yd Roll 48/cs

MSRP: \$118.99 USD

130702 Natural 2.75" x 30 yd Roll 48/cs
MSRP: \$118.99 USD

130705 Big Green 2.75" x 30 yd Roll 48/cs

130707 Big Black 2.75" x 30 yd Roll 48/cs

130709 Big Orange 2.75" x 30 yd Roll 48/cs

130713 Big Maroon 2.75" x 30 yd Roll 48/cs

130715 Big Lime Green 2.75" x 30 yd Roll 48/cs

130721 Big Navy 2.75" x 30 yd Roll 48/cs

Elastic Bandage

High quality cotton and elastic for support and compression.

050101 Beige 2" x 5 yd Roll 10/cs
MSRP: \$36.99 USD

050103 Beige 4" x 5 yd Roll 10/cs

MSRP: \$52.99 USD

050102 Beige 3" x 5 yd Roll 10/cs
MSRP: \$44.99 USD

050104 Beige 6" x 5 yd Roll 10/cs

MSRP: \$67.99 USD

TAPING ACCESSORIES



1. **Tuffner® Quick Drying Spray** Easy to use tape adherent that helps tape and wraps stay in place better and longer
 170211 Tuffner Quick Drying Spray 4 oz aerosol spray 12/cs MSRP: \$8.75 USD 170210 Tuffner Quick Drying Spray 10 oz aerosol spray 12/cs MSRP: \$11.99 USD
2. **M Tac™** A non-aerosol pre-tape spray that helps MTape® and MWrap® stick better
 131301N 8 oz pump spray 12/cs MSRP: \$8.99 USD
3. **Tape & Tuffner® Remover** Removes pre-tape spray and adhesive particles from skin, countertops and taping benches
 200500 4 oz pump citius spray 12/cs MSRP: \$13.99 USD 200115 10 oz aerosol spray 12/cs MSRP: \$10.99 USD
 200103 1 quart 12/cs MSRP: \$37.99 USD 200104 1 gallon MSRP: \$76.99 USD
4. **Heel & Lace Pads** Thin layer of foam to help reduce friction inside shoes, works great with lubricating ointment
 080201 3"x4" pre-cut pads 3/16" thick 1,000 count MSRP: \$150.99 USD
 080202 3"x3" perforated pads 1/16" thick 2,000 count with dispenser MSRP: \$61.99 USD


TAPING SCISSORS

1. **Bandage Scissors** Features an angled blade with blunt ends for quick removal of bandages, tape and more.
 020301 one scissors MSRP: \$9.99 USD
2. **Bandage Shears** Features paramedic-grade stainless-steel construction to cut through straps, padding, clothing and tape.
 020302 one shears MSRP: \$9.99 USD
3. **Emergency Shears** Features pointed high-carbon steel blades sharp enough to cut off face masks.
 020351 one shears MSRP: \$32.99 USD
4. **M Cutter™** Safely slice through stiff athletic tape and wraps.
 100201 one cutter 12/cs MSRP: \$12.99 USD 100202 replacement blades 10/pk MSRP: \$56.99 USD
5. **Economy Tape Cutter** Make removing MTape® or other athletic tape easy.
 200201 one cutter 24/cs MSRP: \$9.99 USD 200301 replacement blades 10/pk MSRP: \$14.99 USD
6. **Kinesiology Scissors** Ideal for precise cutting of kinesiology tape without fraying the edges.
 020328 one scissors MSRP: \$36.99 USD
7. **Super PRO 11 Scissors** Made from a high carbon, hammer forged chromed steel with a small, easy to maneuver safety tip.
 01111 one scissors MSRP: \$68.99 USD
8. **Super PRO 11T Scissors** Made from a high carbon, hammer forged chromed steel. Advanced fluoropolymer non-stick coating.
 01112 one scissors MSRP: \$80.99 USD
9. **Super PRO 21 Scissors** Constructed with extra-sharp diamond honed stainless steel blades and comfort handle.
 01211 one scissors MSRP: \$66.99 USD



10. **Super PRO 21T Scissors** Constructed with extra-sharp diamond honed stainless steel blades. Advanced fluoropolymer non-stick coating.
 01212 one scissors MSRP: \$76.99 USD
11. **Super PRO 21L Scissors** Left handed. Constructed with extra-sharp diamond honed stainless steel blades and comfort handle.
 01213 one scissors MSRP: \$36.99 USD

ATHLETIC TRAINER KITS



BUILT IN SPORT. MADE FOR LIFE.®

UNLIMITED CONFIGURABILITY



Introducing Mueller's extended **premium line** of Athletic Trainer Kits, giving you the power to customize: swap, shift, configure and reconfigure your Kit just the way you need it.

Mueller® has teamed up with one of the industry leaders in first aid bags and is honored to offer a premium line of athletic trainers' bags to meet your every need. Mueller® by MERET™ is built for the professional - professionals that depend on rapidly and proficiently deploying athletic training supplies during crucial moments. Mueller® by Meret™ products will function for a lifetime under the MERET for Life (M4L) Warranty Policy.

Mueller® by Meret™ is built for you, by you.

Visit muellersportsmed.com for more information



1. Medi Kit™ XEL

13020 Black one bag
MSRP: \$699.99 USD

2. Medi Kit™ Omni

13013 Black one bag
13014 Hi Viz Yellow one bag
MSRP: \$539.99 USD

3. Medi Kit™ G2 AT Briefcase

13012 Black one bag
MSRP: \$199.99 USD

4. Medi Kit™ V.E.R.S.A

13018 Black one bag
MSRP: \$369.99 USD

5. Medi Kit™ First In

13015 Black one bag
13016 Hi Viz Yellow one bag
13017 Red one bag
MSRP: \$139.99 USD

6. Medi Kit™ PPE Pro Pack

13005 Black one bag
13006 Hi Viz Yellow one bag
13007 Red one bag
MSRP: \$39.99 USD

7. Medi Kit™ X2 Fill Pak Module

13011 Black one bag
MSRP: \$165.99 USD

8. Medi Kit™ Fill Pak Module

13008 Black one bag
13009 Hi Viz Yellow one bag
13010 Red one bag
MSRP: \$84.99 USD

9. Medi Kit™ MULE

13019 Black one bag
MSRP: \$559.99 USD

MUELLER HERO® KITS

Mueller® Hero® Athletic Trainer kits are customizable with interchangeable accessories. The kits are made with a rugged lightweight nylon exterior with reflective piping for increased visibility in low light. The panels are water-resistant where it is needed the most. A fully-open design lets you unzip bags completely for easy packing and unpacking of your gear. The premium zippers, handles and pockets are built to last; and we back these kits with a one year warranty.

Your possibilities are endless with the customizing options of the Mueller® Hero® Athletic Trainer Kits!



1. Mueller Hero® Utility

includes: M1 Clear Pocket 5 x 10 (1 ea)

19111 one bag
MSRP: \$231.99 USD

2. Mueller Hero® Scout

19922 one bag
MSRP: \$189.99 USD

3. Mueller Hero® Protege

includes: M2 Clear Pocket 10 x 10 (1 ea)

19118 one bag
MSRP: \$209.99 USD
200118 one bag complete with Mueller Refill Kit
MSRP: \$473.99 USD

4. Mueller Hero® Response

includes: M2 Clear Pocket 10 x 10 (1 ea)

19119 one bag
MSRP: \$405.99 USD
200119 one bag complete with Mueller Refill Kit
MSRP: \$517.99 USD

5. Mueller Hero® Perfect

includes: L2 Divider 10 x 10 (2 ea); L2 Divider 5 x 10 (2 ea); M2 Mesh Pocket 10 x 10 (2 ea); Tie Down Strap (4 ea); Rain Fly (1 ea)

19114 one bag
MSRP: \$569.99 USD

6. Mueller Hero® Unsung

includes: L1 Divider 5 x 5 (2 ea); L1 Divider 10 x 5 (2 ea); M2 Mesh Pocket 10 x 10 (1 ea); M1 Clear Folded Pocket (1 ea)

19113 one bag
MSRP: \$352.99 USD

7. Mueller Hero® Epic

includes: L2 Divider 10 x 10 (2 ea); L2 Divider 15 x 10 (2 ea); M2 Mesh Pocket 10 x 10 (1 ea); M2 Clear Pocket 10 x 10 (2 ea); Tie Down Strap (4 ea); Rain Fly (1 ea)

19115 one bag
MSRP: \$675.99 USD

8. Mueller Hero® Titan

includes: L2 Divider 10 x 10 (2 ea); L2 Divider 15 x 10 (2 ea); M2 Mesh Pocket 10 x 10 (1 ea); M2 Clear Pocket 10 x 10 (2 ea); Tie Down Strap (4 ea); Rain Fly (1 ea)

19116 one bag
MSRP: \$750.99 USD



MUELLER HERO® KIT ACCESSORIES

1. Mueller Hero® L1-5 Divider Kit

includes: L1 Divider 5" x 5" (2 ea), L1 Divider 10" x 5" (2 ea)

19551 one divider kit
MSRP: \$14.99 USD

2. Mueller Hero® L2-10 Divider Kit

includes: L2 Divider 10" x 10" (2 ea); L2 Divider 5" x 10" (2 ea)

19552 one divider kit
MSRP: \$48.99 USD

3. Mueller Hero® L2-10 Divider Kit

includes: L2 Divider 10" x 10" (2 ea), L2 Divider 15" x 10" (2 ea)

19553 one divider kit
MSRP: \$57.99 USD

4. Mueller Hero® M1-5 Zipper Pocket Kit

includes: M1 Zipper clear Pocket 10" x 5" (2 ea)

19443 one divider kit
MSRP: \$9.99 USD

5. Mueller Hero® M2-10 Clear Pocket Kit

includes: M2 Clear Pocket 10" x 10" (2 ea)

19773 one pocket kit
MSRP: \$21.99 USD

6. Mueller Hero® M1-5 Clear Folded Pocket Kit

includes: M1 Clear Folded Pocket 10" x 5" (2 ea)

19330 one pocket kit
MSRP: \$22.99 USD

7. Mueller Hero® M2-10 Mesh Pocket Kit

includes: M2 Mesh Pocket 10" x 10" (2 ea)

19772 one pocket kit
MSRP: \$38.99 USD

8. Mueller Hero® M2-10 Folded Clear Bandage Pocket Kit

includes: M2-10 Folded Clear Bandage Pocket 15" x 10" (2 ea)

19220 one pocket kit (19220)
MSRP: \$57.99 USD

9. Mueller Hero® Tie Down Strap Kit

includes: Tie Down Strap (4 ea)

19010 one strap kit
MSRP: \$14.99 USD

10. Mueller Hero® M2-10 Cartridge Kit

includes: M2 Cartridge 10" x 10" (2 ea)

19990 one cartridge kit
MSRP: \$25.99 USD

11. Mueller Hero® M1-10 Cartridge Kit

includes: M1 Cartridge 10" x 5" (2 ea)

19660 one cartridge kit
MSRP: \$19.99 USD



kit contents may vary depending on product availability

Refill Kit

200750 one kit

Kit Includes:

- 12 - Sterile Adhesive Bandages Regular
- 8 - Sterile Adhesive Bandages X-Large
- 12 - Sterile Adhesive Fingertip Bandages
- 1 - 3 oz tube Lubricating Ointment
- 1 - 4 oz Mueller Hand Sanitizer
- 1 - 4 oz Tuffner® Pre-Tape Spray
- 2 - Mueller KOLD® Instant Cold Pack
- 5 - 1/32 oz. Triple Antibiotic Ointment
- 1 - 3" Self Adhering Bandage roll
- 3 - 1 1/2" Athletic Care® Tape roll
- 14 - 3" x 3" Gauze Pads
- 1 - 2 oz Germa™ Spray Plus
- 1 - 3.5 oz Coolant Cold Spray
- 10 - Nasal Plugs
- 1 - Tape Cutter
- 1 - Elastic Bandage roll
- 3 - pair Exam Gloves
- 1 - MWrap® roll
- 1 - Tape Tote™
- 1 - Gauze roll

MSRP: \$89.99 USD



**FITS UP TO
42" WAIST**

overall dimensions: 12" x 4" x 14"

Sling Bag

16017 one bag - sold empty
MSRP: \$81.99 USD

Prop 65 Warning, see p.89 for details



overall dimensions: 12" x 5" x 4.5"

Fanny Pack

200728 one bag - sold empty
MSRP: \$39.99 USD

Prop 65 Warning, see p.89 for details



Tape Roll Holder

16137 one holder - sold empty
MSRP: \$13.99 USD

Prop 65 Warning, see p.89 for details



Tape Tote

200351 one tote - sold empty
MSRP: \$6.99 USD



overall dimensions: 7" x 3" x 5"

First Aid Sport Kit

200738 one complete kit
MSRP: \$29.99 USD

Kit contents may vary depending on product availability

RECOVERY CARE



BUILT IN SPORT. MADE FOR LIFE.®

RECOVERYCARE REVIVE®

Accelerate recovery, reduce pain, and improve performance with new Mueller® REVIVE® Compression Technology. By stimulating the body's natural circulation with precise levels of pressure, REVIVE®'s innovative compression sleeves, Console, and 30% overlapping chambers promote faster recovery before, during, and after workouts, activity, or just standing on your feet. Gear packs consist of console and leg sleeves only. Backed by a 90-day risk-free trial. This adjustable system optimizes your recovery with a personalized range of therapeutic pressure. The durable, portable, and lightweight unit is hassle-free, simple to clean, and powered by 8-hour battery life.



RecoveryCare REVIVE® M4 Gear Pack

The M4 offers 5 unique treatment programs, including a fully customizable program, and is enhanced with Bluetooth technology to track personalized treatments on your mobile device.

M4RG200SM small **M4RG200LG** large
MSRP: \$1,295.00 USD MSRP: \$1,295.00 USD

M4RG200MD medium
MSRP: \$1,295.00 USD



RecoveryCare REVIVE® M2 Gear Pack

At only 2.2 lbs, the M2 is ideal for on-the-go or at home use.

M2RG200SM small **M2RG200LG** large
MSRP: \$789.00 USD MSRP: \$789.00 USD

M2RG200MD medium
MSRP: \$789.00 USD



RecoveryCare REVIVE® Full Leg Boots

Flush out excess toxins in your legs before and after physical activity to accelerate recovery and improve performance.

RG200SM small **RG200LG** large
MSRP: \$450.00 USD MSRP: \$450.00 USD

RG200MD medium
MSRP: \$450.00 USD



RecoveryCare REVIVE® Half Jacket

Targets the upper body area and the core.

RG400
MSRP: \$600.00 USD



RecoveryCare REVIVE® Calf Boots

Provides precise treatment directly to the calf area of the leg. Great for travel!

RG100
MSRP: \$350.00 USD



RecoveryCare REVIVE® Pants Garment

Quickly eliminates excess toxins throughout your entire lower body, accelerating recovery, reducing pain, and improving performance before and after physical activity.

RG300
MSRP: \$799.00 USD



RecoveryCare REVIVE® Core

Targets treatments in hips, glutes, and lower back.

RG600
MSRP: \$625.00 USD



RecoveryCare REVIVE® Arm Sleeve

Provides precise pressure for faster recovery from upper body activity.

RG500
MSRP: \$225.00 USD



RecoveryCare REVIVE® Bag

RG16200
MSRP: \$150.00 USD

RECOVERYTUBS®

Inflatable ice bath, designed for convenience & ease of use! Perfect for teams or athletes on-the-go. The RecoveryTub® also makes an excellent emergency immersion tub for overheating players on the field.

Each unit includes: tub, rechargeable air pump (with adapters), mesh carry bag for easy travel and storage, instruction sheet, and a puncture repair kit.



RecoveryCare RecoveryTub® Solo
 Holds approximately 48 gallons of water
 Empty Weight: 11 lbs
 Fully Inflated dimensions: 59" x 33.4" x 17.7"
 30037
 MSRP: \$349.99 USD



RecoveryCare RecoveryTub® Team
 Holds approximately 106 gallons of water
 Empty Weight: 18 lbs
 Fully Inflated dimensions: 61" Diameter x 27.5" Deep
 30047
 MSRP: \$449.99 USD

INCLUDED ACCESSORIES

Rechargeable Air Pump
 Dimensions: 4" x 3" x 2.3"
 Maximum output pressure: 4 KPA
 Battery capacity: 4000 mAh, 3.6 V, 14.8 W
 USB Charging: 5 V, 2 A



Mesh Carry Bags



THERM-X = HEAT + COLD + CONTRAST + COMPRESSION

Imagine not having to keep refilling, cleaning and drying out your ice based machine. Imagine being able to maintain the specific temperature you want delivered to the patient, for as long as you want. Then consider the new treatments you can apply with heat, compression, and contrast.

This innovative technology allows the user to select Heat, Cold, or Contrast - temperature range of 34° F to 110° F.

With its digital control, the Therm-X circulates a coolant that can maintain a specific temperature for the entire treatment cycle. Compare this to heat packs or ice-based machine that lose temperature once they are applied.



Therm-X AT Machine

The only device that does hot, cold, or contrast therapy with compression at the same time!

TX0002
MSRP: \$6,495.00 USD

THERM-X GARMENTS

Therm-X garments are bilateral so you don't have to buy lefts and rights. They're one size fits all, and have high-end Velcro to last the life of the garment. The inner face uses a ripstop nylon with a waterproof backing to prevent buildup of mold and mildew after cleaning.



Hip Garment

TX0108
MSRP: \$607.00 USD



Back Garment

TX0105
MSRP: \$454.00 USD



Elbow Garment

TX0103
MSRP: \$465.00 USD



Carrying Case

TX01202
MSRP: \$349.00 USD



Ankle Garment

TX0104
MSRP: \$465.00 USD



Knee Garment

TX0102
MSRP: \$482.00 USD



Shoulder Garment

TX0101
MSRP: \$589.00 USD



Garment Extender Strap Kit

TX0109
MSRP: \$144.00 USD

XL Shoulder Garment

TX0110
MSRP: \$682.00 USD



Dual Hose

TX0208
MSRP: \$570.00 USD



Therm-X Coolant 1 Qt

TX0206
MSRP: \$30.00 USD


Replacement Power Supply

2-1030F
MSRP: \$189 USD

Replacement Umbilical Hose- Single


4-1531
MSRP: \$469 USD

COLD/HOT THERAPY



BUILT IN SPORT. MADE FOR LIFE.®

COLD/HOT THERAPY



EZ Relief™ Compression Sleeve / Migraine Cap
Up to 2X longer lasting relief. Thicker gel provides comfort and longer lasting cold temperature to optimize relief and recovery time

6618 Sleeve
MSRP: \$24.99 USD
6613 Migraine Cap
MSRP: \$24.99 USD



Mueller Kold® Instant Cold Pack
Activates instantly for fast response to an injury

030102 6" x 9" 16-Packs (team)
MSRP: \$41.99 USD
330102 6" x 9" 1-Pack 12/cs (retail-clam)
MSRP: \$6.99 USD
35507 6" x 9" 2-Pack (retail-clam)
MSRP: \$9.99 USD

Cold/Hot Pack - Reusable
Soft gel pack can be used for both hot and cold therapy

030104 small 4.75" x 6" 12-Pack (team)
MSRP: \$15.99 USD
030105 regular 6" x 9" 12-Pack (team)
MSRP: \$22.99 USD
030109 therapist 12" x 14" 1-Pack (team)
MSRP: \$20.99 USD
330104 6" x 9" 1-Pack (retail-clam)
MSRP: \$7.99 USD

Beaded Hot/Cold Pack - Reusable
Gel beads allow therapy pack to conform to the contours of your body

32551 Blue 6-Pack (retail)
MSRP: \$59.99 USD
32561 Pink 6-Pack (retail)
MSRP: \$6.99 USD
32557 Blue 1-Pack (retail-clam)
MSRP: \$9.99 USD



Disposable Bag & Dispenser
For quick cold therapy
030801 10" x 18" bag 1,500 Bag-Roll, 3mil thickness (team)
MSRP: \$194.99 USD
030851 dispenser (team)
MSRP: \$35.99 USD

EZ-Wrap
Easy to use, hold cold packs and pads securely in place
050304 EZ-Wrap Roll 4" x 1000' with handle 6/cs (team)
MSRP: \$101.99 USD

Coolant Cold Spray
Cools skin fast and helps relieve swelling
030202EU 9 oz aerosol spray 12/cs (team)
MSRP: \$17.99 USD



Elastic Cold/Hot Wrap

Helps hold cold/hot pack in place

330112 one wrap, one 6" x 9" cold/hot pack (retail-clam)
MSRP: \$17.99 USD



Cold/Hot Therapy Wrap

Adjustable straps help hold pack in place

330121 one wrap, 2 - 4.75" x 6" cold/hot pack for ankle, knee, wrist, forearm (retail-clam)
MSRP: \$21.99 USD
330122 one wrap, 1 - 6" x 9" cold/hot pack for shoulder, thigh, back (retail-clam)
MSRP: \$25.99 USD



Ice Bag Wrap

Easy to adjust to treat multiple areas

230676 one wrap, 9" reusable ice bag (retail-clam)
MSRP: \$27.99 USD



Gel Brace, Cold Therapy

Provides cold therapy while supporting the ankle

4556 OSFM (retail-clam)
MSRP: \$69.99 USD



Mueller Kold Towel

Unique super-evaporating PVA material to help keep you cool

31117 Blue 17" x 12" (retail) MSRP: \$9.99 USD
32117 Blue 17" x 26" (retail)
33117 Pink 17" x 26" (retail)
34117 Kiwi 17" x 26" (retail) MSRP: \$19.99 USD
30291 Blue XL 52" x 32" (retail) MSRP: \$29.99 USD



Fabric Cold/Hot Pack - Reusable

Soft reusable application of cold or heat therapy


36390 12" x 7.8" 1-Pack (retail-clam)
MSRP: \$18.99 USD



Reusable Ice Bag

unique design won't drip or leak

6621A 9" Ice Bag (retail - box)
MSRP: \$18.99 USD



Flexible Cold/Hot Therapy Pads - Reusable

Unique product that can be used hot or cold

030302 one roll w/dispenser 16" x 28" (team)
MSRP: \$108.99 USD

SPORT ACCESSORIES



BUILT IN SPORT. MADE FOR LIFE.®



No Glare® Glare-Reducing Strips

The "original" glare-reducing strips

81407 strips 144/pk (team-bagged)
MSRP: \$35.99 USD

81607 strips 24/pk (retail-blister)
MSRP: \$5.99 USD

81617 pink strips 12/pk (retail-blister)
MSRP: \$5.99 USD



No Glare® Premium Glare-Reducing Strips

Over the nose for added glare protection

82128 strips 12/pk (retail-blister)
MSRP: \$11.99 USD



No Glare® Glare-Reducing Stick

Black paste applied to your cheekbone for glare protection

140403A .65 oz stick 24/cs (team)
MSRP: \$5.99 USD

140413 .15 oz stick 32-ct container (team)
MSRP: \$49.99 USD

440403A .65 oz stick 12/cs (retail-blister)
MSRP: \$5.99 USD

84137 .15 oz stick 12/cs (retail-blister)
MSRP: \$2.99 USD



Foam Rubber, Adhesive Back

Open cell foam rubber

060251 6" x 2 yd roll 1/8" thick
MSRP: \$22.99 USD

060252 6" x 2 yd roll 1/4" thick
MSRP: \$43.99 USD

060271 variety pack: 4 - 6" x 12" 1/8" thick,
4 - 6" x 12" 1/4" thick
MSRP: \$67.99 USD



Felt

Quality felt for extra cushion or added protection

060101 9" x 36" 1/8" thick sheet
MSRP: \$49.99 USD

060151 6" x 2.5 yd roll 1/8" thick, adhesive
MSRP: \$33.99 USD

060152 6" x 2.5 yd roll 1/4" thick, adhesive
MSRP: \$70.99 USD

060161 variety pack: 4 - 6" x 12" 1/8"
thick, adhesive, 4 - 6" x 12" 1/4" thick
adhesive
MSRP: \$57.99 USD




Dermal Pads

Shock absorber for tender skin areas and bruises

040401 4" square 1/8" thick 12/cs
MSRP: \$14.99 USD

040402 4" square 1/2" thick 6/cs
MSRP: \$20.99 USD



More Skin®

Hydrogel that provides moisture to help promote healing

131411 1" square 75/foil bag
MSRP: \$24.99 USD

131412 3" circle 25/foil bag
MSRP: \$74.99 USD

131404 12" x 12" sheet
MSRP: \$22.99 USD



Blister Pads

Thin, flexible, pressure sensitive adhesive foam pads

20257 1.75" x 2.75" foam pads 25/bag
MSRP: \$26.99 USD

200723 1.75" x 2.75" foam pads 5/pack
MSRP: \$11.99 USD



Blister Kit

All the essentials to care for blisters, calluses, abrasions, etc.

includes: 11 foam pressure oval pads, 11 foam pressure donut pads, 1 foam pressure sheet, 6 More Skin® hydrogel pads, 6 small adhesive mesh, 1 large adhesive mesh

200724 one kit
MSRP: \$19.99 USD



Adhesive Mesh

Synthetic fiber with adhesive backing

010901 6" x 10" sheet 15/pack
MSRP: \$68.99 USD



MStrips™

Individually wrapped sterile bandages

230407 regular 1" x 3" 100/box
MSRP: \$9.99 USD

230607 x-large 2" x 4" 50/box
MSRP: \$13.99 USD

260807 fingertip 100/box
MSRP: \$8.99 USD



B Sharp™

Liquid cleaner for cleaning soles of basketball shoes

020101 1 quart bottle 12/cs (team)
MSRP: \$23.99 USD



Germa™ Spray Plus With Lidocaine

Antiseptic spray helps prevent infection in minor cuts, scrapes and burns

33014 2 oz pump spray 24/cs (team)
MSRP: \$5.99 USD



Whizzer® Cleaner & Disinfectant

Cleaner, fungicide, mildewstat, virucide and deodorizer all in one

230201 1 gallon with 1 pump sprayer bottle
MSRP: \$53.99 USD



Instant Hand Sanitizer

Antiseptic gel that kills 99.9% of germs

046004 4 oz bottle 24/cs
MSRP: \$2.99 USD



Sani-Cloth® Plus

Pre-moistened germicidal wipes

MSRP: \$21.99 USD



EZ-Cleans Plus Kit

includes: 1 pair vinyl gloves, 1 pouch Red Z® solidifier 10g, 1 scoop, 1 SaniZidePlus® Germicidal wipe 8" x 8", 1 red bio-hazard waste bag 24" x 24", 1 white disposal bag 30" x 36", 2 twist ties, 1 p.a.w.s® antimicrobial hand wipes 5" x 8", instructions

30733 one kit 12/cs (team)
MSRP: \$25.99 USD



Gloves and Dispenser

Helps protect your hands from messy clean up

070312 medium nitrile gloves 100/box (team)
MSRP: \$52.99 USD

070313 large nitrile gloves 100/box (team)
MSRP: \$52.99 USD

070351 glove dispenser with 6 large nitrile gloves (team)
MSRP: \$9.99 USD



Quart Water Bottle Carrier

Collapsible heavy-duty plastic design for durability

919000 red - holds 6 quart bottles

919111 black - holds 6 quart bottles
MSRP: \$48.99 USD



Quart Water Bottle

Heavyweight, long-lasting BPA free plastic water bottle with easy grip

919129M natural with white straw cap - red logo

919129MB natural with white straw cap - black logo
MSRP: \$6.99 USD

919339M natural with sure shot squeeze cap - red cap & logo

919369MB natural with sure shot squeeze cap - black cap & logo
MSRP: \$7.99 USD




Sterile Gauze Pads

12-ply individually wrapped sterile gauze pads for general first aid

070101 2" x 2" 100/box (team)
MSRP: \$11.99 USD

070102 3" x 3" 100/box (team)
MSRP: \$15.99 USD

070103 4" x 4" 100/box (team)
MSRP: \$21.99 USD



Nasal Sponges

Individually wrapped, fiber free sponges absorb blood to control bleeding

191252 5-Pack (team)
MSRP: \$23.99 USD

191251 25-Pack Jar (team)
MSRP: \$129.99 USD



Cotton Nasal Plugs

Made of soft, non linting, 100% virgin cotton to help control common nosebleeds

191221 300-Pack (team)
MSRP: \$12.99 USD



Tuffner® Quick Drying Spray

Easy to use tape adherent that helps tape and wraps stay in place better and longer

170211 Tuffner Quick Drying Spray 4 oz aerosol spray 12/cs
MSRP: \$8.75 USD

170210 Tuffner Quick Drying Spray 10 oz aerosol spray 12/cs
MSRP: \$11.99 USD



Stickum™ Spray

Spray for hands and equipment to help improve grip

190710 4 oz aerosol spray 12/cs (team)
MSRP: \$8.99 USD



Stickum™ Paste

Designed as a paste to stay on athletes' fingers longer

190801 4 oz jar (team)
MSRP: \$11.99 USD



Stickum™ Grip Powder

Contains rosin in a convenient powder form

490751 1.25 oz shaker (retail)
MSRP: \$12.99 USD



Rosin Blend

Use on hands or soles of your shoes to improve traction

180201A 4 oz shaker (team)
MSRP: \$7.99 USD



Athletic Chalk

Magnesium carbonate in powdered or solid block form

a 2 oz shaker 24/cs (team)
MSRP: \$11.99 USD

070501 8 - 2 oz bars (team)
MSRP: \$33.99 USD

370502 2 oz shaker (retail-blister)
MSRP: \$13.99 USD



Foot Powder

Helps absorb moisture, keeping your feet dry while reducing friction

060301N 4 oz shaker 12/cs (team)
MSRP: \$5.99 USD



Vitamin C Dietary Supplement

Chewable lemon flavored 100 mg tablets

220202 1,000 count jar 6/cs (team)
MSRP: \$80.99 USD



Ball Cleaner

Liquid cleaner used to clean basketballs and footballs

020201 1 quart bottle 12/cs (team)
MSRP: \$12.99 USD



Massage Lotion

Emollient formula provides smooth glide

130801 16 oz bottle 12/cs (team)
MSRP: \$18.99 USD

130803 1 gallon w/pump top (team)
MSRP: \$35.99 USD



Petroleum Jelly

Classic lubricant for friction prone areas

160201 1 lb jar 12/cs (team)
MSRP: \$25.99 USD

160202 5 lb jar 6/cs (team)
MSRP: \$87.99 USD

160203 25 lb drum (team)
MSRP: \$185.99 USD



Lubricating Ointment

Non-friction ointment containing lanolin

120201N 3 oz tube 12/cs (team)
MSRP: \$5.99 USD

120202N 1 lb jar 12/cs (team)
MSRP: \$13.99 USD

120203N 5 lb jar 6/cs (team)
MSRP: \$56.99 USD

120204N 25 lb drum (team)
MSRP: \$225.99 USD



Lube Stick™ For Runners

Provides a barrier between the skin and shoes, pads and equipment

120205N 2.25 oz stick 12/cs (team)
MSRP: \$6.99 USD

420206N .60 oz stick (retail)
MSRP: \$6.99 USD

ATHLETE PROTECTION

ATHLETE PROTECTION



BUILT IN SPORT. MADE FOR LIFE.®

ATHLETE PROTECTION



Diamond Pad Extended Knee
High-impact diamond padding technology for your knee
58206B Black XXXL (team-bagged)
58200-58205 Black XS-XXL (retail-bagged)
(Qty 2, sizing: K1)
MSRP: \$39.99 USD



Diamond Pad Elbow Sleeve
High-impact diamond padding technology for your elbow
75701-75704 Black SM-XL (retail-bagged)
75711-75714 White SM-XL (retail-bagged)
(Qty 1, sizing: E10)
MSRP: \$24.99 USD



Padded Elbow Sleeve
Unique buttress style pad and contour design for flexible protection
417S/M Black SM/MD (retail-clam)
417L/XL Black LG/XL (retail-clam)
417-2X/3X Black XXL-3XL (retail-clam)
(Qty 1, sizing: E11)
MSRP: \$19.99 USD



Strapguard™
Strapguard Mouthguard helps reduce the risk of dental injury during the toughest of sports with a strap to hold secure to your facemask easily
131007 Clear **131013** Gold
131015 Red **131016** Kelly Green
131017 Roayl Blue **131018** Black
131031 Navy Blue **131032** White
131033 Maroon **131034** Purple
131035 Orange **131036** Forest
131037 Pink
MSRP: \$2.99 USD



Muellerguard™
Muellerguard Mouthguard helps reduce the risk of dental injury during the toughest of sports
131005 Clear **131024** Kelly Green
131025 Red **131027** Royal Blue
131028 Black **131029** Gold
131041 Navy Blue **131042** White
131043 Maroon **131044** Purple
131045 Orange **131046** Forest
131047 Pink
MSRP: \$2.99 USD



Youth Mouthguards
Designed to help protect youths' mouth during contact sports
131020 youth with strap
131030 youth without strap
MSRP: \$2.99 USD



Lip & Mouthguard
Lip & Mouthguard helps provide protection to your mouth during contact sports
431019 Black OSFM
MSRP: \$5.99 USD

ATHLETE PROTECTION SIZING CHARTS

MEASURE AROUND KNEE
Standing with a slight bend in your knee, measure the circumference directly around the center of your knee.

K1 - Measure Around Knee

XS	SM	MD	LG	XL	XXL	3XL
10"-12"	12"-14"	14"-16"	16"-18"	18"-20"	20"-22"	22"-24"
25-30 cm	30-35 cm	35-40 cm	40-45 cm	45-50 cm	50-55 cm	55-60 cm

MEASURE AROUND ELBOW
With a slight bend in your elbow, measure the circumference directly around the center of your elbow. If ordering a pair, measure both elbows and find the average.

E10 - Measure Around Elbow

SM	REG	LG	XL
8"-11"	10"-13"	12"-15"	14"-17"
20.3-27.9 cm	25.4-33.0 cm	30.4-38.1 cm	35.5-43.2 cm

E11 - Measure Around Elbow

SM/MD	LG/XL	XXL/XXXL
7"-12"	12"-15"	15"-17"
17-30 cm	30-38 cm	38-43 cm

THE NEXT GENERATION OF PROTECTION

Proprietary Occipital Harness

- Secure fit
- Increased comfort
- Fully adjustable



Medical Grade Polycarbonate

- Crack/shatter proof
- Improved visibility



Anatomically Contoured Design

- Superior fit
- Less obtrusive design



Energy Displacement

- Broadly disperses energy from impact



XPE Foam Padding Kit

- 2 pad thicknesses for:
- Customized fit
- Increased comfort
- Added protection
- Heat management



CUSTOMIZABLE • SHATTER PROOF
INCREASED COMFORT • CONTOURED DESIGN



Face Guard

Helps protect against common impact related sport injuries

81457 1 guard with strap and padding (team-bagged)
MSRP: \$69.99 USD

84457 1 guard with strap and padding (retail-clam)
MSRP: \$69.99 USD

83457 replacement strapping system
MSRP: \$13.99 USD

82457 replacement padding kit
MSRP: \$5.99 USD



TUBULAR COMPRESSION BANDAGES BY PRO ORTHOPEDIC

PRO Orthopedic Tubular Knit Compression Bandages are the high end go-to 'Cut to Fit Compression' bandage option. These bandages offer the best performance value in the industry while out performing more expensive competitors. Made of a tubular, heavy knit material that is very durable. With the high spandex content, it maintains a preferred amount of compression without stretching-out over time. The material has moisture wicking features. It absorbs minimal moisture and is fast drying, providing a cooling effect appreciated by athletes in warm environments. Ideal for securing protective padding, equipment, or ice bags. Bandages are machine washable and can be re-used.



PRO Orthopedic Tubular Knit Compression Bandages
Available in Black or White 11 yard rolls (2", 3", 3.5", 4", 5", 6", 7", and 8")



- 22111 2"x 11 yd roll Black
MSRP: \$47.99 USD
- 23111 3"x 11 yd roll Black
MSRP: \$64.99 USD
- 23511 3.5"x 11 yd roll Black
MSRP: \$51.99 USD
- 24111 4"x 11 yd roll Black
MSRP: \$52.99 USD
- 24511 4.5"x 11 yd roll Black
MSRP: \$55.99 USD
- 25111 5"x 11 yd roll Black
MSRP: \$58.99 USD
- 26111 6"x 11 yd roll Black
MSRP: \$63.99 USD
- 27111 7"x 11 yd roll Black
MSRP: \$74.99 USD
- 28111 8"x 11 yd roll Black
MSRP: \$78.99 USD



- 22191 2"x 11 yd roll White
MSRP: \$47.99 USD
- 23191 3"x 11 yd roll White
MSRP: \$64.99 USD
- 23591 3.5"x 11 yd roll White
MSRP: \$51.99 USD
- 24191 4"x 11 yd roll White
MSRP: \$52.99 USD
- 24591 4.5"x 11 yd roll White
MSRP: \$55.99 USD
- 25191 5"x 11 yd roll White
MSRP: \$58.99 USD
- 26191 6"x 11 yd roll White
MSRP: \$63.99 USD
- 27191 7"x 11 yd roll White
MSRP: \$74.99 USD
- 28191 8"x 11 yd roll White
MSRP: \$78.99 USD

SKY™ ANKLE STABILIZER

HEALTHY ANKLES?

KEEP THEM THAT WAY!

The first of its kind, Mueller has done it again; revolutionizing the game, by preventing sprains and keeping you in top form. Ideal for securing protective padding, equipment, or ice bags. This patented, low-profile ankle stabilizer lowers the risk of a significant inversion injury when landing a vertical leap, as seen commonly in basketball and volleyball players

Velcro® secure straps functionality

Contoured left/right specific design to align with proper anatomy

Curved shape mimics the tape application

Low profile design for minimal bulk in shoe


Semi-rigid stabilizer allows for normal foot flexion and functionality



SKY™ Ankle Stabilizer
MSRP: \$74.99 USD

Left	Right	US Shoe			UK Shoe			EURO Shoe	
		Size	Men	Women	Size	Men	Women	Size	
49631	49641	XXS	3-5	4-6	XXS	2.5-4.5	2-4	XXS	36-38
49632	49642	XS	5-7	6-8	XS	4.5-6.5	4-6	XS	38-40
49633	49643	SM	7-9	8-10	SM	6.5-8.5	6-8	SM	40-42
49634	49644	MD	9-11	10-12	MD	8.5-10.5	8-10	MD	42-44

ANKLE BRACES & SUPPORTS



BUILT IN SPORT. MADE FOR LIFE.®

OUR
BEST *JUST GOT*
BETTER
THE ONE® ANKLE BRACE PREMIUM

**Anatomically Shaped
& Perforated Tongue**
*Better fit with heat &
moisture management*

Low Profile
*Comfortable fit inside
performance footwear*

Figure 8 Strapping
*Mimics traditional taping
techniques*

Grip Treatment
Less slippage means more support

Direct Injected 'SpeedFit' Lacing
MSM Exclusive
*Significantly easier and faster on/off
Engineered eyelet placement enhances fit/function
Greater durability
Flat laces are less binding*

B³ Technology
Breathable Ballistic Boot
MSM Exclusive
*Heat and moisture management
More efficient drying
Conforms to foot shape better*



360° Comfort Band
MSM Exclusive
*Enhanced comfort
Arch support*

Seamless Sole
*Better fit
Greater comfort
More durable*



ADVANCED ANKLE SUPPORTS

Constructed from nylon-vinyl material, these braces are semi-rigid. One of their best features is that you can tighten this brace according to your comfort level. What's more, they're more flexible than rigid braces.

	The One® Ankle Brace Premium	Hg80® Premium Soft Ankle Brace w/Straps	ATF®3 Ankle Brace	Soft Ankle Brace w/ Straps	Adjust-to-Fit Ankle Brace
Sizes Available	XXS-3XL	SM-XXL	XS-XXL	XS-XL	OSFM
Support Level	Maximum	Maximum	Advanced	Advanced	Advanced
Not made with natural rubber latex		×	×	×	
Antimicrobial Treatment*	×	×		×	
Metal Springs					×
Rigid Stays			×		
Lace Up Design	×	×	×	×	×
Adjustable Straps	×	×		×	
Reflective Accents		×			

* The antimicrobial treatment is intended to protect the brace and does not extend protection to the skin.

MAXIMUM



The One® Ankle Brace Premium

- Direct injected speed lace system
- Exclusive B³ body material for superior heat and moisture management

48800B, 48880B-48886B Black XXS-XXXL (team-bagged)

48800, 48880-48884 Black XXS-XL (retail-clam)

(sizing: A1)

MSRP: \$39.99 USD

MAXIMUM



Hg80® Premium Soft Ankle Brace w/Straps

- HydraCinn® high performance fabric
- Antimicrobial treatment

49710-49715 Black XS-XXL (retail-clam)

(sizing: A1)

MSRP: \$47.99 USD

ADVANCED

**ATF@3 Ankle Brace**

- Patented anterior talofibular ligament straps
- Low profile

42374B Black XXL (team-bagged)

42370-42374 Black XS-XXL (retail-clam)

(sizing: A1)

MSRP: \$29.99 USD

ADVANCED

**Soft Ankle Brace w/Straps**

- Low cost, high performance lace up brace
- Strapping system replicates the support of tape

42770-42774 Black XS-XL (team-bagged)

41770-41774 Black XS-XL (retail-clam)

(sizing: A1)

MSRP: \$25.99 USD

ADVANCED

**Adjust-to-Fit Ankle**

- Patented design
- Adjustable side panels for custom fit

4571 Black OSFM (retail-clam)

(sizing: A2)

MSRP: \$24.99 USD

Prop 65 Warning, see p.89 for details

WRAPAROUND ANKLE SUPPORTS

Wraparound ankle supports are a step up from a basic ankle support allowing you to customize the compression for your ankle that is comfortable to you.

	Adjustable Ankle Stabilizer	Adjustable Ankle Support	Mueller Green Adjustable Ankle Support	Ankle Support with Straps	Wraparound Ankle Support
Sizes Available	OSFM	OSFM	OSFM	OSFM	OSFM
Support Level	Advanced	Moderate	Moderate	Moderate	Basic
Not made with natural rubber latex		×	×		×
Antimicrobial Treatment*		×	×		
Rigid Stays	×				
Wraparound Design	×	×	×	×	×
Adjustable Straps	×	×	×	×	

* The antimicrobial treatment is intended to protect the brace and does not extend protection to the skin.



ADVANCED

Adjustable Ankle Stabilizer

- Easy on/off with adjustable compression and support
- Air mesh sole for heat and moisture management

45677 Black OSFM (retail-clam)

(sizing: A5)

MSRP: \$20.99 USD



MODERATE

Adjustable Ankle Support

- Seamless edges for low profile fit
- User adjustable compression

42037 Black OSFM (retail-clam)

42038B Black SM/MD (retail-clam)

(sizing: A4)

MSRP: \$16.99 USD

42039B Black LG/XL (retail-clam)

MODERATE



Mueller Green Adjustable Ankle Support

- Earth friendly recycled body materials (green inside liner)
- Antimicrobial treatment

86511 Black OSFM (retail-box)

(sizing: A3)

MSRP: \$24.99 USD

MODERATE



Ankle Support with Straps

- Fitted slip-on design for snug fit
- User adjustable compression

965SM-965XL Black SM-XL (retail-clam)

(sizing: A1)

MSRP: \$19.99 USD

BASIC



Wraparound Ankle Support

- One size fits most wrap design
- All day comfort

43637 Black OSFM (retail-clam)

(sizing: A2)

MSRP: \$12.99 USD

ANKLE SUPPORTS

These are the simplest type of ankle braces. They provide the least support for activities such as light jogging and walking. However, they are ineffective for intense activities, although they do help in reducing inflammation as the result of injuries.

	OmniForce® Ankle Support A-700	OmniForce® Ankle Support A-100	Hg80® Ankle Support	Elastic Ankle Support	Neoprene Ankle Support
Sizes Available	XS-XXL	XS-XL	XS-XL	SM-XL	SM-XL
Support Level	Moderate	Basic	Moderate	Basic	Moderate
Not made with natural rubber latex	×	×	×		
Antimicrobial Treatment*	×		×		
Gel Pad	×		×		
Step-In Design	×	×	×	×	×
Reflective Accents	×				

* The antimicrobial treatment is intended to protect the brace and does not extend protection to the skin.

MODERATE



OmniForce® Ankle Support A-700

- Advanced 3D flat knit engineering
- Anatomical gel pads for comfort

46611, 42591/46601, 42581 SM 8-10 Women's, 7-9 Men's
 46612, 42592/46602, 42582 MD 10-12 Women's, 9-11 Men's
 46613, 42593/46603, 42583 LG 12-14 Women's, 11-13 Men's
 46614, 42594/46604, 42584 XL 14-16 Women's, 13-15 Men's
 46615/ 46605 XXL 16-18 Women's, 15-17 Men's

(sizing: A1)

MSRP: \$49.99 USD

BASIC



OmniForce® Ankle Support A-100

- High value supportive compression
- Fits left or right

44300B XS 7-8 Women's, 6-7 Men's
 44301 SM 8-10 Women's, 7-9 Men's
 44302 MD 10-12 Women's, 9-11 Men's
 44303 LG 12-14 Women's, 11-13 Men's
 44304 XL 14-16 Women's, 13-15 Men's
 44305B XXL 16-18 Women's, 15-17 Men's

(sizing: A1)

MSRP: \$21.99 USD

MODERATE



4-Way Stretch Ankle Support

- Thermo Reactive Technology
- 360° targeted compression

Fits Left or Right ankle

40008 SM/MD 6-10 Women's, 5-9 Men's

40007 MD/LG 10-14 Women's, 9-13 Men's

40009 LG/XL 14-18 Women's, 13-17 Men's

(sizing: A1)

MSRP: \$15.99 USD

MODERATE



Hg80® Ankle Support

- Lightweight & breathable materials
- Antimicrobial treatment

49910 XS 7-8 Women's, 6-7 Men's

49911 SM 8-10 Women's, 7-9 Men's

49912 MD 10-12 Women's, 9-11 Men's

49913 LG 12-14 Women's, 11-13 Men's

49914 XL 14-16 Women's, 13-15 Men's

(sizing: A1)

MSRP: \$18.99 USD

BASIC



Elastic Ankle Support

- Fitted slip-on design that fits right or left
- Contoured fit

47631 SM 8-10 Women's, 7-9 Men's

47632 MD 10-12 Women's, 9-11 Men's

47633 LG 12-14 Women's, 11-13 Men's

47634 XL 14-16 Women's, 13-15 Men's

(sizing: A1)

MSRP: \$9.99 USD

MODERATE



Neoprene Ankle Support

- Fitted sleeve
- Helps control swelling

41640 XS 7-8 Women's, 6-7 Men's

41641 SM 8-10 Women's, 7-9 Men's

41642 MD 10-12 Women's, 9-11 Men's

41643 LG 12-14 Women's, 11-13 Men's

41644 XL 14-16 Women's, 13-15 Men's

41645B XXL 16-18 Women's, 15-17 Men's

(sizing: A1)

MSRP: \$12.99 USD

RIGID ANKLE BRACES

These braces come in either stirrup or lace-up styles. A stirrup brace, which provides the most protection, lets you move up and down while protecting your ankle from rolling either in or out. These braces are built of a rigid structure that contain plastic panels to protect your ankle from inversion and eversion. A gel or air filled cushion may be part of the brace to help provide comfort to injured ankle.

	Hg80® Premium Hard Shell Ankle Brace	Lite Ankle Brace	Gel Brace, Cold Therapy	Aircast® Ankle Brace
Sizes Available	SM-XL	OSFM	OSFM	OSFM
Support Level	Maximum	Advanced	Maximum	Maximum
Not made with natural rubber latex	×			
Antimicrobial Treatment*	×			
Step-In Design	×	×	×	×
Adjustable Straps	×	×	×	×
Reflective Accents	×			

* The antimicrobial treatment is intended to protect the brace and does not extend protection to the skin.



Hg80® Premium Hard Shell Ankle Brace

- Easy on/off step-in design; no laces
 - HydraCinn® high performance fabric
- 49011-49014 Left SM-XL (retail-clam)
 49021-49024 Right: SM-XL (retail-clam)
 (sizing: A1)
 MSRP: \$54.99 USD



Lite Ankle Brace

- Hinged design aimed at injury prevention or treatment
 - Single strap and easy on/off
- 4552 Black OSFM (retail-clam) 4554 White OSFM (retail-clam)
 (sizing: A2)
 MSRP: \$29.99 USD



Gel Brace, Cold Therapy

- Conforming air chambers
- Removeable gel pads may be frozen for cold therapy

4556 White OSFM (retail-clam)

(sizing: A2)

MSRP: \$69.99 USD



Aircast® Ankle Support

- Pre-inflated air cushions
- Designed for post-injury recovery

4561 White Left OSFM (retail-clam)

4566 White Right OSFM (retail-clam)

(sizing: A2)

MSRP: \$64.99 USD

A1 - According to Shoe Size

	XXS	XS	SM	MD	LG	XL	XXL	XXXL
US MEN	4 - 5	6 - 7	7 - 9	9 - 11	11 - 13	13 - 15	15 - 17	17 - 19
US WOMEN	5 - 6	7 - 8	8 - 10	10 - 12	12 - 14	14 - 16	16 - 18	18 - 20
EURO	36 - 38	38 - 40	40 - 42	42 - 44	44 - 46	46 - 48	48 - 50	50 - 52

A2 - According to Shoe Size

	ONE SIZE
US MEN	7 - 17
US WOMEN	8 - 18
EURO	40 - 50

A3 - According to Shoe Size

	ONE SIZE
US MEN	4.5-12
US WOMEN	6-14
EURO	38 - 40

A4 - According to Shoe Size

	OSFM	SM/MD	LG/XL
US MEN	5 - 13	4 - 7	7 - 13
US WOMEN	6 - 14	5.5 - 8.5	8.5 - 14.5
EURO	38 - 46	37 - 40	40 - 46

A5 - According to Shoe Size

	ONE SIZE
US MEN	3 - 13
US WOMEN	4 - 14
EURO	36 - 46

KNEE BRACES & SUPPORTS




BUILT IN SPORT. MADE FOR LIFE.®

THE ULTIMATE KNIT SUPPORT

OMNIFORCE® KNEE STABILIZER KS-700



Advanced 3D Flat Knit
Creates focused support zones, still allowing the product to stretch and maintain constant compression



Coiled Alloy Tensioners
Provides stabilization on both sides of the knee allowing full range of movement

Breathable Material
Moisture management fabric with superior softness and comfort

Enhanced Length
Better fit and added comfort

Gel Buttress
Provides focused compression to your kneecap

PATENTED

TRIAxIAL HINGE



Single axis and polycentric hinges were designed prior to research that measured the geometric characteristics of knee motion. All other hinges, regardless of misleading names like "multi-centric", are still simple single axis or double axis hinges which are proven to promote a mechanical mismatch that induces brace migration.


NEAR-NORMAL HINGE MOVEMENT

A functional knee brace must produce near-normal joint motion with minimal restrictions. These criteria may be obvious, but clearly many ACL braces fail to meet them.



SHIFT

The posterior shift of the instant center of the knee (I to II). The corresponding femur-tibial contact point.



DISPLACEMENT

Posterior displacement of the sphere centers on the femoral condyles.

COMPARISON OF BRACE HINGE DESIGN

The simplistic movement of single axis and polycentric hinges are incompatible with the complex motion generated by the knee.



SINGLE AXIS HINGE

Hinge design is fixed at the center pivot and exhibits no posterior displacement.



BIAXIAL HINGE

Two centers of rotation allow only minimal posterior displacement.




TRIAxIAL HINGE

The instant center is always posterior of the central axis and reproduces the near-normal posterior pathway of the knee

HINGE BIOMECHANICS

Measurement in the posterior movement of the rotation center in all three types of hinges available are compared to that of a normal knee. Only *Mueller's Triaxial Hinge* properly tracks the joint.

MILLIMETERS OF POSTERIOR MOVEMENT FROM 0° TO 120° OF FLEXION



THINGS TO CONSIDER WHEN CHOOSING A KNEE BRACE

1 WHAT TYPE OF KNEE BRACE DO I NEED?

Knee braces can be categorized into five types: hinges, stabilizers/braces, supports, sleeves and straps. Hinged knee braces contain hinges which allow the knee to move in a way that is close to the knee's natural motion for the most advanced support, while knee stabilizers/braces contain metal springs on both sides of the knee for advanced support. Knee supports can be adjustable and often wrap around the knee; both provide basic support along with warmth and compression, while knee sleeves are sized and can be slipped on over the knee. Basic knee straps can be used to help relieve various knee pains.

2 WHAT STYLE OF KNEE BRACE DO I WANT?

There are two basic styles of knee braces: slip on and wraparound. Slip on braces are slid onto the knee by placing your foot through the brace's opening and pulling up. This type of brace is sized, as the ideal fit is necessary for proper support. Wraparound knee braces can be placed either in front of or behind the knee, depending on the type, and adjusted with the brace's Velcro straps.

3 DO I WANT A KNEE BRACE WITH A CLOSED OR OPEN PATELLA?


Knee braces are made with a closed patella for those who want their knee cap to receive the same compression and support provided to the rest of the knee. Open patella knee braces allow for any pressure on the knee cap to be relieved, and depending on the particular knee brace, may also allow for additional support and/or proper tracking and movement of the knee cap.



PATENTED TRIAXIAL HINGE



DUAL HINGE



SPRINGS

HINGED KNEE BRACES

This type of brace is designed to decrease knee pain and provide stability. They come with hinges that provide sideways stability. A hinged knee brace will keep a knee from bending back too far. With this type of brace, a person is able to determine how much their knee will bend and straighten. Hinged knee braces are able to provide a person with the maximum protection and support.

	Hg80® Premium Hinged Knee Brace	Hg80® Hinged Knee Brace	Pro Level™ Hinged Knee Brace Deluxe	Muellerhinge™ 2100 Knee Brace	Hinged Wraparound Knee Brace	Hinged Knee Brace	Mueller Green Adjustable Hinged Knee Brace
Sizes Available	SM-XXL	SM-XXL	SM-3XL	OSFM	REG-XL	SM-XXL	OSFM
Support Level	Maximum	Maximum	Maximum	Maximum	Maximum	Maximum	Maximum
Patented Triaxial Hinge	✗	✗	✗	✗	✗	✗	
Dual Hinge							✗
Not made with natural rubber latex	✗	✗					✗
Antimicrobial Treatment*	✗	✗			✗		✗
Patella	Open	Closed	Open	Open	Open	Open	Open
Slide On Design	✗	✗				✗	
Wraparound Design			✗	✗	✗		✗
Adjustable Straps	✗	✗	✗	✗	✗	✗	

* The antimicrobial treatment is intended to protect the brace and does not extend protection to the skin.



Hg80® Premium Hinged Knee Brace

- Patented triaxial hinge
- Advanced patterning and assembly methods

59015B Black XXL (team-bagged)

59011-59014 Black SM-XL (retail-clam)

(sizing: K1)

MSRP: \$89.99 USD



Hg80® Hinged Knee Brace

- Patented triaxial hinge
- Tibial Containment System

54015B Black XXL (team-bagged)

54011-54014 Black SM-XL (retail-clam)

(sizing: K1)

MSRP: \$69.99 USD



Pro Level™ Hinged Knee Brace Deluxe

- Patented heavy-duty triaxial hinge
- Removeable patella buttress

53011B SM 12-14 in 53015B XXL 20-22 in
 53012B MD 14-16 in 53016B XXXL 22-24 in
 53013B LG 16-18 in (sizing: K11)
 53014B XL 18-20 in MSRP: \$116.99 USD



Muellerhinge™ 2100 Knee Brace

- Single 12 inch super-duty patented triaxial hinge
- Designed as a protective measure in collision sports

B2100 Black OSFM (team-bagged)
 (sizing: K6)
 MSRP: \$84.99 USD



Hinged Wraparound Knee Brace

- Patented triaxial hinge
- Sectional wrap design

53138 SM/MD 9-14 in
 53137 MD/LG 12-20 in
 53139 LG/XL 17-26 in
 (sizing: K2)
 MSRP: \$52.99 USD



Hinged Knee Brace

- Patented triaxial hinge
- Fitted neoprene sleeve

56320B XS 10-12 in
 56321 SM 12-14 in
 56322 MD 14-16 in
 56323 LG 16-18 in
 56324 XL 18-20 in
 (sizing: K1)
 MSRP: \$64.99 USD



Mueller Green Adjustable Hinged Knee Brace

- Made from earth-friendly recycled materials (green inside liner)
- Dual hinge pockets for sizing

86455 Black OSFM (retail-clam)
 (sizing: K8)
 MSRP: \$49.99 USD

KNEE STABILIZERS & BRACES

Knee Stabilizers & Braces use supportive steel springs on both sides of the knee to help keep it supported but still allowing the joint to move normally. This type of brace is best used for stiff, sore, or floating knees.

	OmniForce® Knee Stabilizer KS-700	Hg80® Premium Knee Brace	Patella Stabilizer Knee Brace	Wraparound Knee Stabilizer	Hg80® Premium Knee Stabilizer	Hg80® Knee Brace
Sizes Available	SM-XXL	SM-3XL	XS-3XL	SM-XXL	SM-XXL	REG-XL
Support Level	Advanced	Advanced	Advanced	Advanced	Advanced	Advanced
Springs	×	×	×	×	×	×
Not made with natural rubber latex	×	×	×	×	×	×
Antimicrobial Treatment*	×	×		×	×	×
Patella	Closed	Open	Open	Open	Open	Closed
Slide On Design	×	×	×		×	×
Wraparound Design		×				
Adjustable Straps		×	×	×		×

	Self Adjusting Knee Stabilizer	Wraparound Knee Brace Deluxe	Elastic Knee Brace	Elastic Knee Stabilizer	Hg80® Knee Support
Sizes Available	OSFM	SM-3XL	SM/MD-LG/XL	SM/MD-LG/XL	XS-XXL
Support Level	Advanced	Advanced	Moderate	Moderate	Advanced
Springs	×	×	×	×	×
Not made with natural rubber latex			×		×
Antimicrobial Treatment*					×
Patella	Open	Open	Open	Open	Closed
Slide On Design				×	×
Wraparound Design	×	×	×		
Adjustable Straps		×	×		

* The antimicrobial treatment is intended to protect the brace and does not extend protection to the skin.



OmniForce® Knee Stabilizer KS-700

- Advanced flat-knit engineering
- 360° gel buttress and two alloy support stays

(sizing: K1)

57771-57775 SM-XXL (team-bagged)

MSRP: \$79.99 USD

53691-53694 SM-XL (retail-clam)

MSRP: \$94.99 USD



Hg80® Premium Knee Brace

- Mueller exclusive HydraCinn® ultra performance fabric with thigh wrap feature
- Four tensioned alloy support stays

59115B XXL (team-bagged)

59111-59114 SM-XL (retail-clam)

(sizing: K1)

MSRP: \$79.99 USD



Patella Stabilizer Knee Brace

- Longer cut fitted sleeve available down to XS
- 'Guitar Pick' shaped patellar support cut-out and four tensioned alloy support stays

55240B XS (team-bagged)

55246B XXXL (team-bagged)

55241-55245 SM-XXL (retail-clam)

(sizing: K1)

MSRP: \$51.99 USD



Wraparound Knee Stabilizer

- 180° inverted support buttress
- Perforated neoprene for heat & moisture management

57638B-57639B S/M-L/XL (team-bagged)

57637 OSFM (retail-clam)

(sizing: K3)

MSRP: \$29.99 USD



Hg80® Premium Knee Stabilizer

- Mueller exclusive HydraCinn® ultra performance fabric
- Fitted sleeve with two tensioned alloy support stays

59215B Black XXL (team-bagged)

59211-59214 Black SM-XL (retail-clam)

(sizing: K1)

MSRP: \$74.99 USD



Hg80® Knee Brace

- Unique Mueller Tibial Containment System with patellar support strap
- HydraCinn® ultra performance fabric and four tensioned alloy support stays

54115B Black XXL (team-bagged)

54111-54114 Black SM-XL (retail-clam)

(sizing: K1)

MSRP: \$59.99 USD



Self Adjusting Knee Stabilizer

- One size fits most wrap design for easy on/off, may be worn over pants
- Criss-Cross elastic straps automatically adjust support above and below knee

56427 Black OSFM (retail-clam)

(sizing: K7)

MSRP: \$34.99 USD



OmniForce® Knee Stabilizer AKS-500

- Hybrid flat knit sectional wrap design
- Triple coiled alloy support stays

(52178, 52179, 52175 retail-clam)

52178 SM/MD 12-16 in, 30-40 cm,

52179 LG/XL 16-20 in, 40-50 cm

MSRP: \$62.99 USD

52175 XXL/XXXL 20-24 in, 50-60 cm

(sizing: K1)

MSRP: \$64.99 USD



Elastic Knee Stabilizer

- Simple yet effective breathable elastic sleeve with form fitting comfort
- Four tensioned alloy support stays

427S/M-427L/XL Black SM/MD-LG/XL (retail-clam)

(sizing: K5)

MSRP: \$16.99 USD



Hg80® Knee Support

- Skin-like fit in a sized sleeve
- Breathable moisture management with antimicrobial treatment and two tensioned alloy support stays

59910-59915 Black XS-XXL (retail-clam)

(sizing: K1)

MSRP: \$29.99 USD

KNEE SUPPORTS

Knee Supports provide support to the knee without interfering with natural range of motion or immobilize the joint. The level of support is different based on the features present in the support. Some are a basic slide on design that provides comfortable compression while others have adjustable straps for a more customized compression.

	4-Way Adjustable Knee Support	Adjustable Knee Support	OmniForce® Knee Support K-700	OmniForce® Knee Support w/Gel K-300	OmniForce® Knee Support K-100	Wraparound Knee Support	Elastic Knee Support
Sizes Available	OSFM SM/MD; LG/XL	OSFM SM/MD; LG/XL	SM-XXL	XS-XL	XS-XL	OSFM SM/MD; LG/XL	SM-XL
Support Level	Moderate	Basic	Moderate	Moderate	Basic	Basic	Basic
Not made with natural rubber latex	✗	✗	✗	✗	✗	✗	
Antimicrobial Treatment*	✗		✗				
Patella	Open	Open	Closed	Closed	Closed	Open	Closed
Slide On Design			✗	✗	✗		✗
Wraparound Design	✗	✗				✗	
Adjustable Straps	✗	✗				✗	

* The antimicrobial treatment is intended to protect the brace and does not extend protection to the skin.



4-Way Adjustable Knee Support

- Easy on/off wrap design may be worn over pants
- 4-Way 'Spider' strap system allows for internal tensioning and direction of support

51238B-51239B Black SM/MD-LG/XL (team-bagged)

51237 Black OSFM (retail-clam)

(sizing: K3)

MSRP: \$28.99 USD



Adjustable Knee Support

- Wrap design with external 'spider' for directional support adjustment
- Engineered not to bunch behind the knee

57228B-57229B Black SM/MD-LG/XL (team-bagged)

57227 Black OSFM (retail-clam)

(sizing: K3)

MSRP: \$20.99 USD



OmniForce® Knee Support K-700

- Advanced flat-knit engineering
 - 360° gel buttress
- 53331-53335 Black SM-XXL (team-bagged)
 52141-52144 Black SM-XL (retail-clam)
 (sizing: K1)
 MSRP: \$79.99 USD



OmniForce® Knee Support with Gel K-300

- Circular knit engineering
 - 360° gel buttress
- 57330B Black XS (team-bagged)
 57331-57334 Black SM-XL (retail-clam)
 (sizing: K1)
 MSRP: \$35.99 USD



OmniForce® Knee Support K-100

- Circular knit engineering
 - Comfortable for all day wear
- 54400B Black XS (team-bagged)
 54401-54404 Black SM-XL (retail-clam)
 (sizing: K1)
 MSRP: \$24.99 USD



Wraparound Knee Support

- Wider straps for better compression adjustment
 - Not made with natural rubber latex
- 53458B-53459B Black SM/MID-LG/XL (team-bagged)
 53457 Black OSFM (retail-clam)
 (sizing: K3)
 MSRP: \$20.99 USD



Elastic Knee Support

- Basic slip-on, low profile support
 - Breathable sized elastic sleeve
- 55251-55254 Black SM-XL (retail-clam)
 (sizing: K1)
 MSRP: \$14.99 USD

KNEE SLEEVES

Knee sleeves are not the same as knee braces. Knee sleeves provide some support to the knee, but do not hinder its natural range of motion or immobilize the joint in any way. This makes knee sleeves unsuitable for providing support after a knee injury or surgery, but rather they are designed to protect the knee from future injury or risk of further damage to an injured area. This protection is particularly important for athletes that put their knees under pressure every day in activities like running, jumping, or weightlifting.

Knee sleeves also add a valuable compression element which help increases blood flow to the area and reduces pain both during and after a workout as well as any post-workout swelling. The neoprene material that knee sleeves are made of make knee sleeves an excellent option to wear post-injury as they provide some support but still allow the knee to strengthen naturally.

	Open Patella Knee Sleeve	Closed Patella Knee Sleeve
Sizes Available	XS-5XL	XS-5XL
Support Level	Moderate	Moderate
Not made with natural rubber latex	✗	✗
Antimicrobial Treatment*		
Patella	Open	Closed
Slide On Design	✗	✗
Wraparound Design		
Adjustable Straps		

* The antimicrobial treatment is intended to protect the brace and does not extend protection to the skin.



Open Patella Knee Sleeve

- Contoured fitted sleeve
- Not made with natural rubber latex

59130B, 59135B-59138B Black XS, XXL-5XL (team-bagged)

59131-59134 Black SM-XL (retail-clam)

(sizing: K1)

MSRP: \$22.99 USD



Closed Patella Knee Sleeve

- Contoured fitted sleeve
- Not made with natural rubber latex

52640B, 52645B-52648B Black XS, XXL-5XL(team-bagged)

52641-52644 Black SM-XL (retail-clam)

(sizing: K1)

MSRP: \$22.99 USD

KNEE STRAPS

Knee straps are designed to compress the tendons surrounding the knee and help redistribute the tension. There is also some evidence that knee straps can be used to keep tendons properly aligned. Knee straps can be worn above or below the knee, depending on the knee injuries plaguing the athlete.

	Max Knee Strap	Advance Patella Strap	Jumper's Knee Strap	ITB Strap
Sizes Available	OSFM SM/MD; LG/XL	OSFM	OSFM	OSFM
Support Level	Maximum	Advanced	Moderate	Moderate
Not made with natural rubber latex	✗	✗	✗	✗
Antimicrobial Treatment*	✗	✗	✗	✗
Patella	Open	N/A	N/A	N/A
Wraparound Design	✗	✗	✗	✗

* The antimicrobial treatment is intended to protect the brace and does not extend protection to the skin.



MAXIMUM

Max Knee Strap

- Upper strap isolates quadriceps misalignment
- Lower strap isolates patellar tendon

59858-59859 Black SM/MD-LG/XL (retail-clam)

59857 Black OSFM (retail-clam)

(sizing: K4)

MSRP: \$29.99 USD



ADVANCED

Advanced Patella Strap



Advanced Patella Strap

- Floating molded EVA buttress for superior function and comfort
- Integrated bi-lateral strap allows for directional adjustment of tension

52898B Black SM/MD (team-bagged)

52899B Black LG/XL (team-bagged)

52897 Black OSFM (retail-box)

(sizing: K3)

MSRP: \$24.99 USD



MODERATE

Jumper's Knee Strap



Jumper's Knee Strap

- Focused pressure on patellar tendon below the knee
- Helps relieve pain from frequent running, jumping and Osgood-Schlatter's

51997 Red OSFM (retail-clam)

52997 Black OSFM (retail-clam)

53997 Blue OSFM (retail-clam)

54997 Gold OSFM (retail-clam)

55997 Pink OSFM (retail-clam)

(sizing: K9)

MSRP: \$14.99 USD



MODERATE

ITB Strap



ITB Strap


- Helps relieve lateral pain above the knee
- Textured neoprene provides anti-slip effect

57927 Black OSFM (retail-clam)


(sizing: K9)

MSRP: \$22.99 USD


KNEE SIZING CHARTS



MEASURE AROUND KNEE
Standing with a slight bend in your knee, measure the circumference directly around the center of your knee.



MEASURE AROUND 3 IN ABOVE KNEE
Standing with a slight bend in your knee, measure 3 inches above your knee cap and measure the circumference around your leg at that point.



MEASURE AROUND THIGH AT WIDEST POINT
Measure the circumference of the thickest part of your thigh.

MEASURE AROUND CALF AT WIDEST POINT
Measure the circumference of the thickest part of your calf. If ordering a pair, measure both calves and find the average.

K1 - Measure Around Knee

XS	SM	MD	LG	XL	XXL	3XL	4XL	5XL
10"-12"	12"-14"	14"-16"	16"-18"	18"-20"	20"-22"	22"-24"	24"-26"	26"-28"
25-30 cm	30-35 cm	35-40 cm	40-45 cm	45-50 cm	50-55 cm	55-60 cm	60-66 cm	66-71 cm

K2 - Measure Around Knee

REG	LG	XL
12"-16"	16"-20"	20"-24"
30-40 cm	40-50 cm	50-60 cm

K3 - Measure Around Knee

ONE SIZE	SM/MD	LG/XL
12"-20"	9"-15"	17"-26"
30-50 cm	23-38 cm	43-66 cm

K4 - Measure Around Knee

ONE SIZE	SM/MD	LG/XL
12"-18"	9"-12"	12"-18"
30-45 cm	22.8-30.4 cm	30.4-45.7 cm

K5 - Measure Around Knee

SM/MD	LG/XL
12"-16"	16"-20"
30-40 cm	40-50 cm

K6 - Measure Around Knee

ONE SIZE
10"-26"
25-66 cm

K7 - Measure Around Knee

ONE SIZE
14"-20"
35-50 cm

K8 - Measure Around 3 in Above Knee

ONE SIZE
13"-21"
33-53 cm

K9 - Measure Around Knee

ONE SIZE
10"-22"
25-55 cm

K10 - Measure Around 3 in Above Knee

SM	MD	LG	XL	XXL	3XL
10"-13"	13"-15"	15"-17"	17"-19"	20"-21.5"	21.5"-23"
25-33 cm	33-38 cm	38-43 cm	43-48 cm	50-54 cm	54-58 cm

K11 - Measure Around Thigh & Calf

	SM	MD	LG	XL	XXL	XXXL
THIGH	17" - 18" 43 - 45 cm	19" - 20" 48 - 50 cm	21" - 22" 53 - 55 cm	23" - 24" 58 - 60 cm	25" - 26" 63 - 66 cm	27" - 28" 68 - 71 cm
CALF	13" - 14" 33 - 35 cm	14" - 15" 35 - 38 cm	15" - 16" 38 - 40 cm	16" - 17" 40 - 43 cm	17" - 18" 43 - 45 cm	18" - 19" 45 - 48 cm

MUELLER GREEN

EARTH FRIENDLY

BRACES & SUPPORTS

Mueller Green is created from a material made without the use of petroleum, and the green inner lining is manufactured with 67% recycled plastic PET bottles. The spandex binding is made from 73% recycled materials. Even the packaging is recyclable and printed with soy-based inks, bringing our commitment to the earth full circle.

SOFT, DURABLE OUTER FABRIC

CORE MADE WITH
NON-PETROLEUM BASED MATERIALS

LINING MADE OF 67% RECYCLED PET BOTTLES

FROM BOTTLES TO FABRIC: THE GREEN PROCESS

Common PET bottles are washed, ground and extruded into fabric.



1. PET Bottle



2. Flakes



3. Chips



4. Yarn



5. Fabric

Adjustable Hinged
Knee Brace p.49



Adjustable Ankle
Support p.41



Fitted Wrist Brace p.70



Adjustable Back & Abdominal
Support p.76



ELBOW BRACES & SUPPORTS



BUILT IN SPORT. MADE FOR LIFE.®

MUELLER
OMNIFORCE®
 SUPPORTS & STABILIZERS



OmniForce®
 Knee Stabilizer
 KS-700 p.50

OmniForce® Supports & Stabilizers are made of a custom lightweight knit design that helps to support your joints for both performance and everyday activities. The comfortable materials and contoured fit make the Omniforce® line great for extended use, no matter the activity.

OmniForce® 100 Series



OmniForce® Ankle
 Support A-100 p.40



OmniForce® Elbow
 Support E-100 p.61



OmniForce® Knee
 Support K-100 p.53

Reflective threading pattern
 Circular knit provides for a seamless design
 Targeted compression vertically and horizontally

OmniForce® 300 Series



OmniForce® Knee
 Support K-300 p.53

Reflective threading pattern
 Circular knit provides for a seamless design
 Contoured shape helps prevent slipping and bunching
 Targeted compression vertically and horizontally

OmniForce® 700 Series



OmniForce® Ankle
 Support A-700 p.40



OmniForce® Elbow
 Support E-700 p.61



OmniForce® Knee
 Support K-700 p.53

Flat knit creates focused support zones while still allowing the product to stretch and maintain constant compression
 Contoured shape helps prevent slipping and bunching
 Targeted compression vertically and horizontally
 Gel Pads for added comfort

ADJUSTABLE ELBOW SUPPORTS

The adjustable straps on these elbow supports provide a customized fit while allowing full range of motion.

	Adjustable Elbow Support	Elbow Support
Sizes Available	OSFM	OSFM
Support Level	Advanced	Moderate
Not made with natural rubber latex		
Antimicrobial Treatment*	✗	
Slide On Design	✗	
Wraparound Design		✗
Adjustable Straps	✗	✗

* The antimicrobial treatment is intended to protect the brace and does not extend protection to the skin.



Adjustable Elbow Support

- Slip on design for easy adjustment
- HydraMesh® panel for heat and moisture management

75217 Black OSFM (retail-clam)

(sizing: E1)

MSRP: \$20.99 USD



Elbow Support

- Adjustable straps
- Comfortable, soothing compression

4521 Black OSFM (retail-clam)

(sizing: E2)

MSRP: \$14.99 USD



4-Way Stretch Elbow Support

- Thermo Reactive Technology
- 360° targeted compression

Fits Left or Right Elbow

Measure circumference of the elbow

(retail-clam)

70058 SM/MD 7-10 in

70057 MD/LG 10-13 in

(sizing: E2)

MSRP: \$14.99 USD

ELBOW SUPPORTS & SLEEVES

	OmniForce® Elbow Support E-700	OmniForce® Elbow Support E-100	Hg80® Elbow Support	Elastic Elbow Support	Professional Elbow Sleeve	Elbow Sleeve	Turf Sleeve	Performance Sleeve
Sizes Available	SM-XXL	SM-XXL	SM-XXL	SM-XL	YTH, REG, LG	SM-XL	LG, XL	OSFM
Support Level	Moderate	Basic	Moderate	Basic	Moderate	Basic	Basic	Basic
Not made with natural rubber latex	×	×	×					×
Antimicrobial Treatment*	×		×					
Slide On Design	×	×	×	×	×	×	×	×

* The antimicrobial treatment is intended to protect the brace and does not extend protection to the skin.



OmniForce® Elbow Support E-700

- Advanced 3-D flat knit construction
- Anti-Microbial treatment

78881-78885 Black SM-XXL (team-bagged)

72361-72364 Black SM-XL (retail-clam)

(sizing: E5)

MSRP: \$46.99 USD



OmniForce® Elbow Support E-100

- Circular knit
- Sized fit

72205B Black XXL (team-bagged)

72201-72204 Black SM-XL (retail-clam)

(sizing: E5)

MSRP: \$21.99 USD

MODERATE



Hg80® Elbow Support

- Super thin/lightweight
- Gel pads provide targeted compression

79911-79915 Black/Silver SM-XXL (retail-clam)
 (sizing: E5)
 MSRP: \$18.99 USD

BASIC



Elastic Elbow Support

- Lightweight elastic knit
- Custom top weave helps prevent rolling

74181-74184 Black SM-XL (retail-clam)
 (sizing: E7)
 MSRP: \$11.99 USD

MODERATE



Professional Reversible Elbow Sleeve

- Compression, support, and protection
- Popular among pro football players

410Y Black Youth (retail-clam)
 410REG Black Regular (retail-clam)
 410LG Black Large (retail-clam)
 (sizing: E6)
 MSRP: \$19.99 USD

BASIC



Elbow Sleeve

- Contoured fitting sleeve
- Available in 4 sizes

414SM-414XL Black SM-XL (retail-clam)
 (sizing: E7)
 MSRP: \$15.99 USD

BASIC



Turf Sleeve

- Abrasion protection for forearm and elbow
- Sold as a pair

4111LG-4111XL LG-XL (retail-clam)
 (sizing: E7)
 MSRP: \$11.99 USD

BASIC



Performance Sleeve

- Seamless design; no stitching lengthwise
- Engineered cuffs help prevent bunching

70007 Black OSFM 70017 White OSFM
 70027 Pink OSFM 70037 Red OSFM
 70047 Blue OSFM
 (sizing: E8)
 MSRP: \$17.99 USD

TENNIS ELBOW SUPPORTS

	Hg80® Premium Tennis Elbow	Tennis Elbow w/Support Gel Pad	Tennis Elbow Support
Sizes Available	SM/MD; LG/XL	OSFM	OSFM
Support Level	Advanced	Advanced	Moderate
Not made with natural rubber latex	✗	✗	
Antimicrobial Treatment*	✗	✗	
Slide On Design	✗	✗	✗
Adjustable Straps	✗	✗	✗
Padding	✗	✗	

* The antimicrobial treatment is intended to protect the brace and does not extend protection to the skin.



Hg80® Premium Tennis Elbow Brace

- Highly effective 'Counter-Force' engineering
- Mueller exclusive HydraCinn® Fabric

79018-79019 Black SM/MD-LG/XL (retail-clam)

(sizing: E9)

MSRP: \$24.99 USD

Tennis Elbow Support w/Gel Pad

- Gel pad targeted compression
- Anti-Microbial treatment

70207 Black OSFM (retail-clam)

(sizing: E4)

MSRP: \$15.99 USD



Tennis Elbow Support

- Wide strap
- Durable neoprene blend


819 Black OSFM (retail-clam)

(sizing: E1)


MSRP: \$13.99 USD

820 Beige OSFM (retail-clam)

ELBOW SIZING CHARTS



MEASURE AROUND ELBOW
With a slight bend in your elbow, measure the circumference directly around the center of your elbow. If ordering a pair, measure both elbows and find the average.



MEASURE AROUND BICEP
With a slight bend in your elbow, slightly flex and measure the circumference of the thickest part of your bicep. If ordering a pair, measure both biceps and find the average.

E1 - Measure Around Elbow

ONE SIZE

7"-14"

17-35 cm

E2 - Measure Around Elbow

ONE SIZE

10"-18"

25-45 cm

E3 - Measure Around Elbow

ONE SIZE

8"-16"

20-40 cm

E4 - Measure Around Elbow

ONE SIZE

10"-14"

25-35 cm

E5 - Measure Around Elbow

SM MD LG XL XXL

7"-9" 9"-11" 11"-13" 13"-15" 15"-17"

18-23 cm 23-28 cm 28-33 cm 33-38 cm 38-43 cm

E6 - Measure Around Elbow

YTH REG LG

9"-11" 11"-14" 14"-17"

22-27 cm 27-35 cm 35-43 cm

E7 - Measure Around Elbow

SM MD LG XL

7"-10" 10"-12" 12"-13" 13"-15"

17-25 cm 25-30 cm 30-33 cm 33-38 cm

E8 - Measure Circumference of Bicep

ONE SIZE

14"-18"

35.6-45.7 cm

E9 - Measure Around Forearm

SM/MD LG/XL

7"-11" 11"-15"

17-27 cm 27-38 cm

MUELLER HG80® PREMIUM BRACES & SUPPORTS

Hg80® Premium braces and supports maximize comfort and support. Made of a Mueller-exclusive HydraCinn® fabric, you can expect high-performance, moisture-wicking material with an antimicrobial barrier.*

UP TO 30% LIGHTER THAN BASIC HG80 PRODUCTS

REFLECTIVE ACCENTS FOR INCREASED VISIBILITY
IN LOW-LIGHT CONDITIONS

SOFT, LIGHTWEIGHT MATERIALS

IMPROVED STRETCH OF MATERIALS

IMPROVED SHAPE TO BETTER FIT NATURAL
SHAPE OF YOUR BODY



Hg80® Premium
Knee Brace p.50



Hg80® Premium
Knee Stabilizer p.51



Hg80® Premium
Soft Ankle Brace
w/Straps p.36



Hg80® Premium
Tennis Elbow Brace p.63



Hg80® Premium
Wrist Brace p.67



Hg80® Premium
Hinged Knee
Brace p.48

WRIST BRACES & SUPPORTS



BUILT IN SPORT. MADE FOR LIFE.®

WRIST BRACES & SUPPORTS

Wrist Braces & Stabilizers are a more rigid brace to help create stability to your wrist or hand. These feature stabilizers such as metal springs, rigid stays or metal spoons to help minimize mobility while promoting healing of your wrist or hand.

	Hg80® Premium Wrist Brace	Mueller Green Fitted Wrist Brace	Wrist Brace w/Splint	Wrist Stabilizer	Thumb Stabilizer	Thumb Stabilizer	Wrist Brace
Sizes Available	SM/MD; LG/XL	SM/MD; LG/XL	OSFM	SM/MD; LG/XL	OSFM	OSFM	OSFM
Support Level	Advanced	Maximum	Maximum	Maximum	Maximum	Advanced	Advanced
Not made with natural rubber latex	✗	✗		✗	✗		
Antimicrobial Treatment*	✗	✗					
Steel Springs	✗						✗
Rigid Stays		✗		✗	✗		
Metal Spoon			✗	✗			
Wraparound Design			✗				✗
Slide On Design	✗	✗		✗	✗	✗	
Adjustable Straps	✗	✗	✗	✗	✗	✗	✗

* The antimicrobial treatment is intended to protect the brace and does not extend protection to the skin.

ADVANCED



Hg80® Premium Wrist Brace

- Mueller exclusive HydraCinn® material for advanced heat and moisture management
- Compact yet supportive

Fits Left or Right Elbow

79718-79719 Black SM/MD-LG/XL (retail-clam)

(sizing: W11)

MSRP: \$37.99 USD

MAXIMUM



Mueller Green Fitted Wrist Brace

- Made from recycled materials (green inside liner)
- Contoured fit is great for slender wrists

86271 Black Right SM/MD (retail-box)

86272 Black Left SM/MD (retail-box)

(sizing: W2)

MSRP: \$24.99 USD

86273 Black Right LG/XL (retail-box)

86274 Black Left LG/XL (retail-box)

MAXIMUM



Wrist Brace w/Splint

- Removeable splint
- One size fits most; fits either wrist

300 Black OSFM (retail-clam)
 (sizing: W9)
 MSRP: \$21.99 USD

MAXIMUM



Wrist Stabilizer

- Maximum support; extra-long 6.5" removeable alloy splint
- Fits left or right wrist; two sizes

307-308 Black SM/MD-LG/XL (retail-clam)
 (sizing: W2)
 MSRP: \$25.99 USD

MAXIMUM



Thumb Stabilizer

- Breathable mesh fabric
- Anti-microbial treatment

42717 Black OSFM (retail-clam)
 (sizing: W7)
 MSRP: \$24.99 USD

ADVANCED



Thumb Stabilizer

- Shock absorbing foam
- Metal splint for stability

(sizing: W7)
 MSRP: \$19.99 USD

ADVANCED



Wrist Brace

- Coiled alloy stays provide support
- One size wrap design

222 Black OSFM (retail-clam)
 (sizing: W10)
 MSRP: \$19.99 USD

Prop 65 Warning, see p.89 for details

WRIST SUPPORTS

	Neoprene Wrist Sleeve	Wrist Support Wrap	Elastic Wrist Support w/Loop	Hg80® Wrist Support	Elastic Wrist Support	Compression Gloves
Sizes Available	SM-XL	OSFM; SM/MD; LG/XL	OSFM	SM-XL	REG; LG	SM-LG
Support Level	Moderate	Moderate	Basic	Moderate	Basic	Moderate
Not made with natural rubber latex		✗		✗		
Antimicrobial Treatment*		✗		✗		
Wraparound Design	✗	✗	✗			
Slide On Design				✗	✗	✗
Adjustable Straps	✗	✗	✗			✗

* The antimicrobial treatment is intended to protect the brace and does not extend protection to the skin.



Neoprene Wrist Sleeve

- Sized neoprene sleeve provides cushioned compression
- Fits either wrist

400SM-400XL Black SM-XL (retail-clam)

(sizing: W5)

MSRP: \$17.99 USD



Wrist Support Wrap

- Simple, quick and easy to adjust
- Perforated neoprene provides breathable comfort

75138 Black SM/MD (retail-clam)

75139 Black LG/XL (retail-clam)

4505 Black OSFM (retail-clam)

(sizing: W3)

MSRP: \$12.99 USD



Elastic Wrist Support w/Loop

- Easy to apply with one hand
- Basic support fits either wrist

961 Black OSFM (retail-clam)

(sizing: W8)

MSRP: \$8.99 USD



Hg80® Wrist Support

- Comfortable compressive support in four sizes; fits either wrist
- Antimicrobial treatment

79111-79114 Black SM-XL (retail-clam)

(sizing: W4)

MSRP: \$18.99 USD



Elastic Wrist Support

- Contoured breathable elastic knit
- Reinforced cuffs for durability

76058 Black Regular (retail-clam)

76059 Black Large (retail-clam)

(sizing: W6)

MSRP: \$9.99 USD



Compression Gloves (pair)

- Comfortable, breathable elastic fabric
- Fingerless design

6901 Beige SM (retail-box)


6903 Beige MD (retail-box)

6904 Beige LG (retail-box)

(sizing: W4)

MSRP: \$16.99 USD

WRISTS SIZING CHARTS



MEASURE AROUND WRIST
With your wrist straight, measure the circumference directly around your wrist.

W1 - Measure Around Wrist

SM/MD	LG/XL
6"-8"	8"-9"
15-20 cm	20-22 cm

W2 - Measure Around Wrist

SM/MD	LG/XL
5"-8"	8"-10"
12-20 cm	20-25 cm

W3 - Measure Around Wrist

OSFM	SM/MD	LG/XL
6"-12"	6"-8"	6"-12"
15-30 cm	15-20 cm	15-30 cm

W4 - Measure Around Wrist

SM	MD	LG	XL
6.5"-7.5"	7.5"-8.5"	8.5"-9.5"	9.5"-10.5"
16-19 cm	19-21 cm	21-24 cm	24-26 cm

W5 - Measure Around Wrist

SM	MD	LG	XL
6"-7"	7"-8"	8"-9"	9"-10"
15-17 cm	17-20 cm	20-22 cm	22-25 cm

W6 - Measure Around Wrist

REG	LG
5"-6.5"	6.5"-8"
12-16.5 cm	16.5-20 cm

W7 - Measure Around Wrist

ONE SIZE
5.5"-10.5"
13-26 cm

W8 - Measure Around Wrist

ONE SIZE
6"-12"
15-30 cm

W9 - Measure Around Wrist

ONE SIZE
6"-9"
15-22 cm

W10 - Measure Around Wrist

ONE SIZE
7"-10"
17-25 cm

W11 - Measure Around Wrist

SM/MD	LG/XL
6"-8"	8"-9"
15-20 cm	20-22 cm

BACK BRACES & SUPPORTS

BACK BRACES & SUPPORTS



BUILT IN SPORT. MADE FOR LIFE.®

BACK BRACES & SUPPORTS

	Lumbar Back Brace w/ Removable Pad	Adjustable Back Brace	Back Support w/Suspenders	Multi Purpose Back Brace	Adjustable Back & Abdominal Support	Mueller Green Back & Abdominal Support	Waist Support
Sizes Available	REG; Plus Size	OSFM	OSFM	SM-XXL	OSFM	OSFM	OSFM
Support Level	Maximum	Advanced	Advanced	Advanced	Maximum	Maximum	Basic
Not made with natural rubber latex				x	x	x	x
Antimicrobial Treatment*				x		x	
Steel Springs	x	x		x			
Stays			x				
Padding	x						
Adjustable Tension	x	x	x	x	x	x	
Tapered Cut	x	x					x

* The antimicrobial treatment is intended to protect the brace and does not extend protection to the skin.



Lumbar Back Brace w/Removable Pad

- Optional insert allows for focused support over the spine
- Multi-layer engineering to prevent rolling or bunching

255 Black Regular (retail-clam)

MSRP: \$42.99 USD

Prop 65 Warning, see p.89 for details

64179 Black Plus Size (retail-clam)

MSRP: \$59.99 USD

(sizing: B2)



Adjustable Back Brace

- 8 coiled alloy lumbar support stays
- Easy adjust contoured fit

4581 Black OSFM (retail-clam)

(sizing: B3)

MSRP: \$28.99 USD



Back Support w/Suspenders

- Adjustable/removeable suspenders
- Breathable material with elastic tensioning bands

252 Black OSFM (retail-clam)
 (sizing: B4)
 MSRP: \$35.99 USD



Adjustable Back & Abdominal Support

- Unique one-handed cinch strap for effective adjustments
- Helps promote proper posture

68147 Black OSFM (retail-clam)
 (sizing: B5)
 MSRP: \$34.99 USD
 Prop 65 Warning, see p.89 for details



Mueller Green Adjustable Back & Abdominal

- Unique one-handed cinch strap for effective adjustments
- Made from recycled materials (green inside liner)

86741 Black OSFM (retail-box)
 (sizing: B5)
 MSRP: \$49.99 USD




Waist Support

- Conforming support to help slim the waist
- Single large band for easy on and off

68127 Black OSFM (retail-clam)
 (sizing: B3)
 MSRP: \$28.99 USD

BACK SIZING CHARTS



MEASURE AROUND WAIST
 Standing in a relaxed position, measure around your waist at the narrowest point.

B1 - Measure Around Waist

SM	MD	LG	XL	XXL
28"-32"	32"-36"	36"-40"	40"-44"	44"-48"
71-81 cm	81-91 cm	91-101 cm	101-111 cm	111-121 cm

B2 - Measure Around Waist

REG	PLUS SIZE
28"-50"	50"-70"
71-127 cm	127-178 cm

B3 - Measure Around Waist

ONE SIZE
28"-50"
71-127 cm

B4 - Measure Around Waist

ONE SIZE
26"-44"
66-111 cm

B5 - Measure Around Waist

ONE SIZE
32"-51"
81.2-129.5 cm

GRADUATED COMPRESSION

BUILT IN SPORT. MADE FOR LIFE.®

SMART PERFORMANCE & RECOVERY



BENEFITS

- ▶ Graduated compression helps improve blood & oxygen circulation
- ▶ Helps reduce muscle vibration, which may cause soreness and fatigue

FEATURES

- ▶ Antimicrobial* garment odor control
- ▶ Moisture management material
- ▶ Extremely durable fabric
- ▶ Breathable performance
- ▶ Mueller exclusive HYDRASHADE™ UV protection fabric UPF 50+

GRADUATED COMPRESSION

Graduated compression helps improve blood and oxygen circulation while helping to reduce muscle vibration which could cause soreness and fatigue. Equipped with Mueller Exclusive HYDRASHADE™ UV Protection fabric UPF 50+

Compression Levels

Performance 20-30 mmHg

Recovery 15-20 mmHg



Graduated Compression Socks

Performance

20-30 mmHg

(retail-box)

- 42020 XS 5-7 Women's, 4-6 Men's
- 42021 SM 7-9 Women's, 6-8 Men's
- 42022 MD 9-11 Women's, 8-10 Men's
- 42023 LG 11-13 Women's, 10-12 Men's
- 42024 XL 13-15 Women's, 12-14 Men's
- 42025 XXL 15-17 Women's, 14-16 Men's

(sizing: G1)

MSRP: \$49.99 USD

Graduated Compression Socks

Recovery

15-20 mmHg

(retail-box)

- 43020 XS 5-7 Women's, 4-6 Men's
- 43021 SM 7-9 Women's, 6-8 Men's
- 43022 MD 9-11 Women's, 8-10 Men's
- 43023 LG 11-13 Women's, 10-12 Men's
- 43024 XL 13-15 Women's, 12-14 Men's
- 43025 XXL 15-17 Women's, 14-16 Men's

(sizing: G1)

MSRP: \$49.99 USD

Graduated Compression Ankle Socks

20-30 mmHg

47020-47025 Black

48020-48025 White

(retail-box)

- 47020, 48020 XS 5-7 Women's, 4-6 Men's
- 47021, 48021 SM 7-9 Women's, 6-8 Men's
- 47022, 48022 MD 9-11 Women's, 8-10 Men's
- 47023, 48023 LG 11-13 Women's, 10-12 Men's
- 47024, 48024 XL 13-15 Women's, 12-14 Men's
- 47025, 48025 XXL 15-17 Women's, 14-16 Men's

(sizing: G1)

MSRP: \$23.99 USD




Graduated Compression Calf Sleeves Performance

45020-45025 Black XS-XXL pair (retail-box)
(sizing: G3)
MSRP: \$44.99 USD



Graduated Compression Arm Sleeves Performance



72020-72025 Black XS-XXL pair (retail-box)
(sizing: G4)
MSRP: \$39.99 USD



Graduated Compression Leg Sleeves Performance

46020-46026 Black XS-XXL pair (retail-box)
(sizing: G2)
MSRP: \$69.99 USD

GRADUATED COMPRESSION SIZING CHARTS



G1 - According to Shoe Size

	XS	SM	MD	LG	XL	XXL
US MEN	4 - 6	6 - 8	8 - 10	10 - 12	12 - 14	14 - 16
US WOMEN	5 - 7	7 - 9	9 - 11	11 - 13	13 - 15	15 - 17

G3 - Measure Calf Circumference at Widest Point

XS	SM	MD	LG	XL	XXL
9"-11"	11"-13"	13"-15"	15"-17"	17"-19"	19"-21"
22-27 cm	27-33 cm	33-38 cm	38-43 cm	43-48 cm	48-53 cm

G2 - Measure Thigh Circumference at Widest Point

XS	SM	MD	LG	XL	XXL	XXXL
15"-19"	17"-21"	19"-23"	21"-25"	23"-27"	25"-29"	27"-31"
38-48 cm	43-53 cm	48-58 cm	53-63 cm	58-68 cm	63-73 cm	68-78 cm

G4 - Measure Bicep Circumference at Widest Point

XS	SM	MD	LG	XL	XXL
8"-10"	10"-12"	12"-14"	14"-16"	16"-18"	18"-20"
20-25 cm	25-30 cm	30-35 cm	35-40 cm	40-45 cm	45-50 cm

G2 - Measure Calf Circumference at Widest Point

XS	SM	MD	LG	XL	XXL	XXXL
13"-14"	14"-15"	15"-16"	16"-17"	17"-18"	18"-19"	19"-20"
33-35 cm	35-38 cm	38-40 cm	40-43 cm	43-45 cm	45-48 cm	48-50 cm

SPECIALTY BRACES & SUPPORTS



BUILT IN SPORT. MADE FOR LIFE.®

SPECIALTY BRACES & SUPPORTS



Shoulder Support

- Multi-strap adjustment system provides a stable custom fit while allowing full range of motion
- Recommended for muscle strains, ligament sprains, injuries to rotator cuff, impingement syndrome, arthritis and bursitis of the shoulder

69530 Black XS (team-bagged)

69538 Black SM/MD (team-bagged)

69539 Black LG/XL (team-bagged)

69535 Black XXL (team-bagged)

(sizing: S1)

MSRP: \$103.99 USD



Arm Sling

- Thumb loop to help hold in place
- Adjustable strap

4601 Blue OSFM (retail-clam)

MSRP: \$14.99 USD

Prop 65 Warning, see p.89 for details



Thigh Support

- Wraparound design for custom fit
- Soft neoprene blend retains body heat to help improve circulation

4491 Black OSFM (retail-clam)

(sizing: S2)

MSRP: \$21.99 USD



Thigh Sleeve

- Fitted, uniform compression
- Soft neoprene blend retains body heat to help improve circulation

444SM-444XL Black SM-XL (retail-clam)

(sizing: S3)

MSRP: \$21.99 USD

MODERATE



Multi-Directional Calf Wrap

- Customizable wrap
- Multi-directional fastening provides uniform down pressure on shin

Measure around calf
(retail-clam)

40018 SM/MD 11-16 in

40019 LG/XL 16-20 in

40015B XXL/XXXL 18-24 in

(sizing: S4)

MSRP: \$29.99 USD

OmniForce® Plantar Fascia Sock

- Provides targeted compression to the arch and heel
- All-day comfort and relief

40137 Black OSFM (retail-clam)

40138 Black SM/MD (retail-clam)

40139 Black LG/XL (retail-clam)

(sizing: S7)

MSRP: \$20.99 USD



Arch Support

- Wraparound design
- Gel pad for added comfort and targeted pressure

46027 Black OSFM (retail-clam)

(sizing: S5)

MSRP: \$14.99 USD



Pro Heel Cups

For high-impact sports, helps divert shock away from the heel

971REGA Green Regular pair (retail-clam)

971LGA Green Large pair (retail-clam)

(sizing: S6)

MSRP: \$13.99 USD

SPECIALTY SIZING CHARTS

S1 - Measure Around Center of Chest

XS	SM/MD	LG/XL	XXL
26" - 34"	32" - 40"	38" - 46"	44" - 52"
66-86 cm	81-101 cm	96-116 cm	106-132 cm

S2 - Measure Around Thigh

ONE SIZE
15" - 35"
38-88 cm

S3 - Measure Around Thigh

SM	MD	LG	XL
21" - 23"	23" - 25"	25" - 27"	27" - 29"
53-58 cm	58-63 cm	63-68 cm	68-73 cm

S4 - Measure Around Calf

REG	LG
11" - 16"	16" - 20"
27-40 cm	40-50 cm

S5 - Measure Around Arch

ONE SIZE
up to 18"
up to 45 cm

S6 - By Body Weight


REG	LG
under 175 lbs	over 175 lbs
under 79.4 kilos	over 79.4 kilos

S7 - According to Shoe Size

	OSFM	SM/MD	LG/XL
US MEN	6 - 12	3 - 7	11 - 16
US WOMEN	7 - 13	4 - 8	12 - 17
EURO	39 - 45	36 - 40	44 - 49

QUENCH[®] GUM

QUENCH[®] GUM



BUILT IN SPORT. MADE FOR LIFE.[®]

QUENCH® GUM

Need that little bit of 'pick-me-up' during sports? Quench® Gum is the answer. It's the fat-free thirst quenching gum for active people. Available in four taste-tingling flavors: Lemon, Fruit, Orange and Grape. Quench is a mouthwatering explosion!



Quench® Gum Variety Box Shelf Talker
Display box with variety boxes
17200 12 boxes per display
MSRP: \$0.15/piece \$0.25/2-piece USD



Tub-O-Quench
300 assorted pieces of gum per tub
17100 4 tubs per case
MSRP: \$0.15/piece \$0.25/2-piece USD



Quench® Gum Variety Bag
2.4 oz variety bag of assorted gum flavors

17600 12 variety bags per case
17700 36 variety bags with three clip strips per case
MSRP: \$0.15/piece \$0.25/2-piece USD



Quench® Gum Team Sports Bottle
Premium quart size sports bottle with assorted gum flavors
1919329QE 6 bottles per case
MSRP: \$9.99 bottle



Quench® Gum Sports Team Bucket
200 assorted pieces of gum per bucket
170576 4 buckets per case
MSRP: \$0.15/piece \$19.99/bucket



Quench® Gum 10-Stick Packs
Mouthwatering double flavors
170143-T Orange Fruit Punch 12/Tray
170153-T Strawberry Watermelon 12/Tray
170157-T Double Raspberry 12/Tray
MSRP: \$1.59/pack USD

170145-T Orange Fruit Punch 24/Tray
170155-T Strawberry Watermelon 24/Tray
170159-T Double Raspberry 24/Tray



Quench® Gum 10-Stick Pack Shelf Talker
Countertop double flavor display: 48-10 stick packs per shelf talker
170144 Orange Fruit Punch
170147 Double Raspberry
MSRP: \$1.59/pack USD
170146 Strawberry Watermelon

4-Way Adjustable Knee Support.....	52	F	Fabric Cold/Hot Pack, Reusable.....	24
A			Face Guard.....	31
Adhesive Mesh.....	26		Fanny Pack.....	17
Adjustable Ankle Stabilizer.....	38		Felt.....	26
Adjustable Ankle Support.....	38		First Aid Sport Kit.....	17
Adjustable Back & Abdominal Support.....	73		Flexible Cold/Hot Therapy Pads, Reusable.....	24
Adjustable Back Brace.....	72		FlexLight™ Tape.....	7
Adjustable Elbow Support.....	60		Foam Rubber, Adhesive Back.....	26
Adjustable Knee Support.....	52		Foot Powder.....	28
Adjust-to-Fit Ankle.....	37	G		
Advanced Patella Strap.....	55		Gel Brace, Cold Therapy.....	24, 43
Aircast® Ankle Support.....	43		Germa™ Spray Plus With Lidocaine.....	27
Ankle Support with Straps.....	39		Gloves and Dispenser.....	27
Arch Support.....	80		Graduated Compression Ankle Socks.....	76
Arm Sling.....	79		Graduated Compression Arm Sleeves.....	77
ATF®3 Ankle Brace.....	37		Graduated Compression Calf Sleeves.....	77
Athletic Care® Athletic Trainers' Tape.....	6		Graduated Compression Leg Sleeves.....	77
Athletic Care® Porous Tape.....	6		Graduated Compression Socks.....	76
Athletic Chalk.....	28	H		
Athletic Trainer Kit, Refil Kit.....	17		Heel Cups.....	80
B			Heel & Lace Pads.....	12
Back Support w/Suspenders.....	73		Hg80® Ankle Support.....	41
Ball Cleaner.....	28		Hg80® Elbow Support.....	62
Bandage Scissors.....	12		Hg80® Hinged Knee Brace.....	48
Bandage Shears.....	12		Hg80® Knee Brace.....	51
Beaded Hot/Cold Pack, Reusable.....	23		Hg80® Knee Support.....	51
Blister Kit.....	26		Hg80® Premium Braces & Supports.....	65
Blister Pads.....	25		Hg80® Premium Hard Shell Ankle Brace.....	42
B Sharp™.....	27		Hg80® Premium Hinged Knee Brace.....	48
C			Hg80® Premium Knee Brace.....	50
Closed Patella Knee Sleeve.....	54		Hg80® Premium Knee Stabilizer.....	51
Cold/Hot Pack, Reusable.....	23		Hg80® Premium Soft Ankle Brace w/Straps.....	36
Cold/Hot Therapy Wrap.....	24		Hg80® Premium Tennis Elbow Brace.....	63
Compression Gloves.....	69		Hg80® Premium Wrist Brace.....	67
Coolant Cold Spray.....	23		Hg80® Wrist Support.....	69
Cotton Nasal Plugs.....	28		Hinged Knee Brace.....	49
D			Hinged Wraparound Knee Brace.....	49
Dermal Pads.....	26	I		
Diamond Pad Elbow Sleeve.....	30		Ice Bag, Reusable.....	24
Diamond Pad Extended Knee.....	30		Ice Bag Wrap.....	24
Disposable Bag & Dispenser.....	23		Instant Hand Sanitizer.....	27
E			ITB Strap.....	55
Economy Tape Cutter.....	12	J		
Elastic Ankle Support.....	41		Jumper's Knee Strap.....	55
Elastic Bandage.....	11	K		
Elastic Cold/Hot Wrap.....	24		Kinesiology Scissors.....	12
Elastic Elbow Support.....	62		Kinesiology Tape Continuous Rolls.....	9
Elastic Knee Stabilizer.....	51	L		
Elastic Knee Support.....	53		Lip & Mouthguard.....	30
Elastic Wrist Support.....	69		Lite Ankle Brace.....	42
Elastic Wrist Support w/Loop.....	69		Lube Stick™ For Runners.....	28
Elbow Sleeve.....	62		Lubricating Ointment.....	28
Elbow Support.....	60		Lumbar Back Brace w/Removable Pad.....	72
Emergency Shears.....	12	M		
EZ-Cleans Plus Kit.....	27		Massage Lotion.....	28
EZ-Wrap.....	23		Max Knee Strap.....	55
EZ-Relief Compression Sleeve.....	23		M Cutter™.....	12
EZ-Relief Migraine Cap.....	23			

Medi Kit™ Fill Pack Module.....	14
Medi Kit™ MULE.....	14
Medi Kit™ Omni.....	14
Medi Kit™ V.E.R.S.A.....	14
Medi Kit™ XEL.....	14
Meret™ Kit Accessories.....	14
MLastic® Tape.....	8
More Skin®.....	26
MStrips™.....	26
M Tac™.....	12
MTape®.....	6
Mueller Green Adjustable Ankle Support.....	39
Mueller Green Adjustable Back & Abdominal.....	73
Mueller Green Adjustable Hinged Knee Brace.....	49
Mueller Green Fitted Wrist Brace.....	67
Muellerguard.....	30
Mueller Hero® Epic.....	15
Mueller Hero® Kit Accessories.....	16
Mueller Hero® Kits.....	16
Mueller Hero® Perfect.....	15
Mueller Hero® Protege.....	15
Mueller Hero® Response.....	15
Mueller Hero® Scout.....	15
Mueller Hero® Titan.....	15
Mueller Hero® Unsung.....	15
Mueller Hero® Utility.....	15
Muellerhinge™ 2100 Knee Brace.....	49
Mueller Kold® Instant Cold Pack.....	23
Mueller Kold Towel.....	24
MWrap®.....	11
N	
Nasal Sponges.....	28
Neoprene Ankle Support.....	41
Neoprene Wrist Sleeve.....	69
No Glare® Glare-Reducing Stick.....	26
No Glare® Glare-Reducing Strips.....	26
No Glare® Premium Glare-Reducing Strips.....	26
O	
OmniForce® Ankle Support A-100.....	40
OmniForce® Ankle Support A-700.....	40
OmniForce® Elbow Support E-100.....	61
OmniForce® Elbow Support E-700.....	61
OmniForce® Knee Stabilizer KS-700.....	50
OmniForce® Knee Support K-100.....	53
OmniForce® Knee Support K-700.....	53
OmniForce® Knee Support with Gel K-300.....	53
OmniForce® Plantar Fascia Sock.....	80
Open Patella Knee Sleeve.....	54
P	
Padded Elbow Sleeve.....	30
Patella Stabilizer Knee Brace.....	50
Performance Sleeve.....	62
Perform Plus™ Tape.....	6
Perform™ Tape.....	6
Petroleum Jelly.....	28
Professional Elbow Sleeve.....	62

Pro Level™ Hinged Knee Brace Deluxe.....	49
ProStrips®.....	10
PTape®.....	10
Q	
Quart Water Bottle.....	27
Quart Water Bottle Carrier.....	26
Quench® Gum.....	81
R	
Recoil™ Tape.....	8
Rosin Blend.....	28
S	
Sani-Cloth® Plus.....	27
Self Adjusting Knee Stabilizer.....	51
Shoulder Support.....	79
Sling Bag.....	17
Soft Ankle Brace w/Straps.....	37
Sterile Gauze Pads.....	28
Stickum™ Grip Powder.....	28
Stickum™ Paste.....	28
Stickum™ Spray.....	28
Strapguard.....	30
Stretch MTape® Premium.....	7
T	
Tape Roll Holder.....	17
Tape Tote.....	17
Tape & Tuffner® Remover.....	12
Tapewrap® Premium.....	8
Tear-Light Tape®.....	7
Tennis Elbow Support.....	63
Tennis Elbow Support w/Gel Pad.....	63
The One® Ankle Brace Premium.....	35
Thigh Sleeve.....	79
Thigh Support.....	79
Thor™ Tape.....	8
Thumb Stabilizer.....	68
Tuffner® Quick-Drying Spray.....	11
Turf Sleeve.....	62
V	
Vitamin C Dietary Supplement.....	28
W	
Waist Support.....	73
Whizzer® Cleaner & Disinfectant.....	27
Wraparound Ankle Support.....	39
Wraparound Knee Stabilizer.....	51
Wraparound Knee Support.....	53
Wrist Brace.....	68
Wrist Brace w/Splint.....	68
Wrist Stabilizer.....	68
Wrist Support Wrap.....	69

Some products in this catalog contain the following warnings:

 **WARNING**
Cancer and Reproductive Harm
www.P65Warnings.ca.gov

DGER SPORTIN



VISIT US AT:
MUELLERSPORTSMED.COM

IG GOODS CO., INC.



For questions or assistance, please call 800-356-9522 to connect with our Mueller Family.



Mueller Sports Medicine, Inc. One Quench Drive, PO Box 99, Prairie du Sac, WI 53578 USA • 1-800-356-9522 • 608-643-8530
ORDER FAX 1-800-852-4334 • FAX 608-643-2568 • sportcare@muellersportsmed.com
muellersportsmed.com

999014E-6.23 ©2023 Mueller Sports Medicine, Inc.